



VELTEK ASSOCIATES, INC.

TECHNICAL DATA FILES



HYPO-CHLOR[®] Neutral 0.25% & 0.52%

Neutralized Sodium Hypochlorite 0.25% & 0.52% Formulated with Water for Injection
Sterile Pharmaceutical Cleanroom Formula

Claims, product registrations, and regulatory requirements may vary based on local, regional, and/or national laws and regulations. Please contact your local sales representative for information specific to your region.

15JAN2025

Copyright © 2025 Veltek Associates, Inc.

DECON, EPA

All Rights Reserved. Confidential and Proprietary Property of Veltek Associates, Inc. No part of this document may be reproduced or disclosed, in whole or in part, without the specific written permission of Veltek Associates, Inc.



TABLE OF CONTENTS

PRODUCT OVERVIEW 3

 Quality and Manufacturing 3

 Features and Benefits 3

 Product Uses 4

 ABCD Cleanroom Introduction System® 5

SIMPLEMIX® 5

 The SimpleMix Technology Alternative 5

 The Activator 5

 SimpleMix Directions for Use 6

CHEMICAL SPECIFICATIONS 10

ORDERING INFORMATION 11

LABELING 12

 VAI’s Product Label Colors 12

 Product Labeling 13

TESTING SUMMARY 22

 HYPO-CHLOR Neutral Efficacy Summary 22

ADDITIONAL DOCUMENTATION 27



HYPO-CHLOR Neutral Family of Products

PRODUCT OVERVIEW

HYPO-CHLOR Neutral products are an effective, neutral pH sodium hypochlorite solution formulated with Water for Injection (WFI) at 0.25% and 0.52% concentrations. EPA registered **HYPO-CHLOR Neutral** products are available in our convenient, ready-to-use, SimpleMix® System, designed for all pharmaceutical, biopharmaceutical, biotechnology, health care, medical device, and diagnostic manufacturing cleaning and disinfecting rotations that demand a neutral pH sodium hypochlorite solution adequate for maintaining a clean and critical environment. This neutralized solution can be used as an improved and enhanced sodium hypochlorite and hard surface Disinfectant, Fungicide, Germicide, Virucide*, and Sporicide on numerous cleanroom surfaces. By lowering the pH of the sodium hypochlorite, the concentration of the hypochlorous acid increases and surface corrosivity decreases.

HYPO-CHLOR Neutral products are manufactured via aseptic fill at 0.2 microns into gamma irradiated sterile components in an ISO 5 cleanroom (Grade A). Each lot of **HYPO-CHLOR Neutral** is sterility tested according to current USP Compendium and is completely traceable. **HYPO-CHLOR Neutral** products are delivered each time with lot specific analytical and sterility data, tested to current USP compendium.

HYPO-CHLOR Neutral 0.25% and 0.52% concentrations are available sterile in a 16 fl. oz. trigger spray, a 1 gallon container, and a 200 L drum. All container sizes come in our one-step, ready-to-use, SimpleMix System that allows for exact and fresh solutions each and every time without handling the activator or sodium hypochlorite solution. Each sterile container is individually double bagged and packaged in two liner bags using the ABCD Cleanroom Introduction System®.

When used per SimpleMix System directions, **HYPO-CHLOR Neutral** products are neutralized via the attached activator when the product is required for use.

**When used as directed*

Quality and Manufacturing

- Formulated with Water for Injection
- Filled in an ISO 5 cleanroom (Grade A)
- Filtered at 0.2 microns
- Components are air washed with 0.2 micron filtered air before assembly
- Aseptically filled into sterile components via gamma irradiation
- Lot sterility tested according to current USP compendium
- Completely traceable from start to finish

Features and Benefits

- EPA Registered
 - EPA Registration Number: 68959-9 and 68959-11
- Each sterile container is double bagged packaged in easy tear bags
- Quadruple bagged in the ABCD Cleanroom Introduction System®

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
15JAN2025 / DECON, EPA

Features and Benefits (Continued)

- Individually labeled with lot number and expiration
- Delivered with lot specific analytical and sterility data, tested to current USP compendium
- Specifically formulated as a sterile cleanroom pharmaceutical formula
- Available in two convenient, one-step, ready-to-use SimpleMix concentrations: 0.25% and 0.52%
- Available in 16 fl. oz., 1 gallon, or 200 L containers
- Comes in a convenient 16 fl. oz. trigger spray that has the option of spray or stream
- Effective for up to 7 days post activation
- Enhanced cleaning applications over a standard sodium hypochlorite solution
- Hard Surface Disinfectant, Fungicide, Germicide, Virucide*, and Sporicide
- Increased cleaning surface compatibility
- Neutralized sodium hypochlorite will significantly reduce corroding, rusting, and pitting of cleanroom surfaces
- Designed for all washable non-porous environmental surfaces

**When used as directed*


Product Uses

- Cleanrooms and controlled areas in health care institutions, biopharmaceutical, pharmaceutical, medical device and diagnostic manufacturing facilities
- Hard non-porous, inanimate surfaces
- Aseptic filling and gowning rooms, general manufacturing areas, laboratories
- Machinery, tools, tables, counters, laminar-flow benches, floors, walls, carts, shelves
- Plastic, glass, vinyl, glazed porcelain, laminates, glazed tiles, stainless steel

This neutral pH sodium hypochlorite solution is compatible with most non-porous hard surface materials while reducing corroding, rusting, and pitting. **HYPO-CHLOR Neutral** 0.25% and 0.52% are safe to use daily on stainless steel surfaces. In addition, testing has demonstrated that **HYPO-CHLOR Neutral** products are effective in deactivating Endotoxins during specialized manufacturing processes.

ABCD Cleanroom Introduction System®

The ABCD Cleanroom Introduction System is a packaging system that allows operators/users to take the package through each level of classified areas by simply removing one bag at a time. Each bag acts as barrier protecting the finished product from becoming a carrier of viable and non-viable contamination. This prevents the need to decontaminate each outer bag prior to entering a cleaner area. In this packaging system, sterilized groups of containers are contained in two outer bags and after each are removed individual containers are each additionally contained in two easy tear bags.



SIMPLEMIX®

The SimpleMix Technology Alternative

Veltek Associates, Inc. has developed the patented SimpleMix System Technology to eliminate measuring and additional containers. It provides for the transfer of the sterile concentrate and sterile water in a sealed container to the aseptic area. The system container is double bag packaged for easy transfer and eliminates all internal and external sterility concerns. The patented SimpleMix System provides a sealed multi-chamber container that when activated mixes the solution to the correct use dilution. The small chamber holds the concentrate, and the bottom reservoir contains the VAI WFI Quality Water. Just open the small chamber cap, push the plunger container completely down until the bottom pops open and the bellows are compressed. The solution and water mix together. 200L SimpleMix systems are activated through a hose and valve system connecting the cubicontainer of concentrate to the VAI® WFI Quality Water. The SimpleMix System design permits the easy transfer of the product to the aseptic manufacturing area without concern for the transfer of contamination.

The Activator

HYPO-CHLOR Neutral 0.25% and 0.52% are activated via the SimpleMix System by neutralizing the pH. The activator is housed in the small chamber suspended above the sodium hypochlorite solution in the large chamber below. Once the SimpleMix plunger is pushed, the activator is mixed into the solution following the SimpleMix directions. Testing has shown that **HYPO-CHLOR Neutral** products are effective as neutralized cleaners for up to 7 days post activation.



SimpleMix Bottle with Activator

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
 Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
 15JAN2025 / DECON, EPA

SimpleMix® Directions for Use

How to use the SimpleMix Trigger Spray

1. To prepare use solution, open cap.
2. Peel off inner seal by grasping tab at far edge and pulling off.
3. Firmly push small, inner container completely down.
4. Replace cap and tighten.
5. Slowly swirl for 15 seconds.
6. Move spray nozzle to open position.
7. Follow directions for use on label.

For pictorial directions, see page 7.

How to use the SimpleMix Gallon Size

1. To prepare use solution, open cap.
2. Peel off inner seal by grasping far edge and pulling off.
3. Firmly push small, inner container completely down.
4. Replace cap and tighten.
5. Slowly swirl for 15 seconds.
6. Open small side spout and peel off inner seal, as above.
7. Pour solution from small side spout onto surfaces to be treated or alternate containers.
8. Follow directions for use on label.

For pictorial directions, see page 8.

How to use the SimpleMix 200 Liter Drum

1. Close all valves.
2. Uncoil hoses.
3. Connect center hose to pump between X and Y.
4. Open valve 1, then valve 2, then valve 4.
5. START pump to empty cubic container.
6. When cubic container is empty, turn pump OFF.
7. Close valve 1 and valve 2.
8. Open valve 6 and valve 5.
9. Re-start pump and mix 15 minutes
10. STOP pump.
11. Close valve 4.
12. To dispense – Open valves 3 and 7. Run pump only when dispensing.
13. Follow directions for use on label.

For pictorial directions, see page 9.



SimpleMix Trigger Spray

SIMPLEMIX® 16 oz/ 473mL Aseptic Mixing System

For the Exact Formulation of 16 oz/ 473 mL Disinfectants and Sporicides

Ready-to-Use Mixing Instructions

- 1) To prepare use solution, open cap.
2) Peel off inner seal by grasping tab at far edge and pulling off.



- 3) Firmly push small, inner container all the way down.



- 4) Replace cap and tighten.



- 5) Slowly swirl for 15 seconds.



- 6) Move spray nozzle to open position.



- 7) Follow directions for use on label.



Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234

Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com

15JAN2025 / DECON, EPA

SimpleMix Gallon

SIMPLEMIX® 1 Gallon/ 3.79L Aseptic Mixing System For the Exact Formulation of 1 Gallon/ 3.79L Size Disinfectants and Sporicides **Ready-to-Use Mixing Instructions**

- 1) To prepare use solution, open cap.
- 2) Peel off inner seal by grasping tab at far edge and pulling off.



4) Replace cap and tighten.



6) Open small side spout and peel off inner seal, as above.



- 3) Firmly push small, inner container all the way down.



5) Slowly swirl for 15 seconds.



7) Pour solution from small side spout onto surfaces to be treated or alternate containers.



8) Follow directions for use on label.

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
15JAN2025 / DECON, EPA

SimpleMix Drum

SIMPLEMIX® 200 L Aseptic Mixing System For Large Scale Aseptic Manufacturing Environments **Ready-to-Use Mixing Instructions**

Remove drum from double-bag packaging.



Remove cubic container from top of drum. 1) Close all valves. 2) Uncoil hoses.



3) Connect center hose to pump between X and Y.



4) Open valve 1, then valve 2, then valve 4.



5) START pump to empty cubic container. 6) When cubic container is empty, turn pump OFF.



7) Close valve 1 and valve 2.



8) Open valve 6 and valve 5.



9) Re-start pump and mix 15 minutes. 10) Stop pump.



11) Close valve 4. 12) To dispense- Open valves 3 and 7. Run pump only when dispensing.



13) Follow directions for use on label.

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
15JAN2025 / DECON, EPA

CHEMICAL SPECIFICATIONS


| HYPO-CHLOR Neutral 0.25% - 0.25% Sodium Hypochlorite formulated with Water for Injection | |
|---|-----------------------|
| Certificate of Analysis | Specifications |
| Assay | 0.25 - 0.31% w/w |
| Litmus paper turns blue | Pass |
| Addition of HCL gives off CL ₂ gas | Pass |
| Yellow flame test | Pass |
| Expiration period | 2 years |

| HYPO-CHLOR Neutral 0.52% - 0.52% Sodium Hypochlorite formulated with Water for Injection | |
|---|-----------------------|
| Certificate of Analysis | Specifications |
| Assay | 0.43 – 0.63% w/w |
| Litmus paper turns blue | Pass |
| Addition of HCL gives off CL ₂ gas | Pass |
| Yellow flame test | Pass |
| Expiration period | 2 years |

| Attached Activator for 0.25% and 0.52% – Neutral Aqueous Solution | |
|--|-----------------------|
| Certificate of Analysis | Specifications |
| Appearance | Colorless |
| pH | 6.5 – 7.5 |
| Specific Gravity | 1.2 – 1.4 |
| Expiration period once activated | 7 days |

ORDERING INFORMATION

| HYPO-CHLOR Neutral 0.25% and 0.52% – Neutralized Sodium Hypochlorite 0.25% and 0.52% Formulated with Water for Injection | | |
|---|--|---------------|
| Part Number | Description | Qty/cs |
| SHC-NPH-0.25-16Z | HYPO-CHLOR Neutral 0.25%, 16 oz SimpleMix, Attached Activator, Attached Trigger, Sterile | 12 |
| SHC-NPH-0.25-02 | HYPO-CHLOR Neutral 0.25%, 1 Gallon SimpleMix, Attached Activator, Sterile | 4 |
| SHC-NPH-0.25-200L | HYPO-CHLOR Neutral 0.25%, 200 L SimpleMix Drum, Attached Activator, Sterile | 1 |
| SHC-NPH-0.52-16Z | HYPO-CHLOR Neutral 0.52%, 16 oz SimpleMix, Attached Activator, Attached Trigger, Sterile | 12 |
| SHC-NPH-0.52-02 | HYPO-CHLOR Neutral 0.52%, 1 Gallon SimpleMix, Attached Activator, Sterile | 4 |
| SHC-NPH-0.52-200L | HYPO-CHLOR Neutral 0.52%, 200 L SimpleMix Drum, Attached Activator, Sterile | 1 |



SHC-NPH-0.52-02



SHC-NPH-0.25-16Z

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
 Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
 15JAN2025 / DECON, EPA



LABELING

VAI's Product Label Colors

| Product Name | Bottle/Can Color | Label Background Color | Bar & User Info Color | Text Color |
|--|------------------------|------------------------|-----------------------|------------|
| DECON-AHOL WFI FORMULA 70% AEROSOL | COOL GREY | PRINTED CAN COOL GREY | | |
| DECON-AHOL WFI FORMULA 70% TRIGGER SPRAY, 1 & 5 GALLON | WHITE | COOL GREY | | |
| DECON-AHOL WFI FORMULA 70% SQUEEZE BOTTLE | WHITE SEMI-TRANSPARENT | COOL GREY | | |
| DECON-AHOL WFI FORMULA 70% ASEPTI-CLEANSE BOTTLE | WHITE SEMI-TRANSPARENT | COOL GREY | | |
| DECON-AHOL WFI FORMULA 60% | WHITE | COOL GREY | | |
| DECON-AHOL WFI FORMULA 91% | WHITE | COOL GREY | | |
| DECON-AHOL FORMULA 99% | WHITE | COOL GREY | | |
| STER-AHOL WFI AEROSOL | WHITE | PRINTED CAN WHITE | | |
| STER-AHOL WFI TRIGGER SPRAY, 1 & 5 GALLON | WHITE | WHITE | | |
| DECON-HAND STERILE | WHITE SEMI-TRANSPARENT | PRINTED BOTTLE | | |
| DECON-HAND NON-STERILE | CLEAR | PRINTED BOTTLE | | |
| DECON-HAND ASEPTI-CLEANSE BOTTLE | WHITE SEMI-TRANSPARENT | WHITE | | |
| STERI-OIL | WHITE | WHITE | | |
| STERI-BUFFER | CLEAR | WHITE | | |
| DECON-PHENE | WHITE | WHITE | | |
| DECON-CYCLE | WHITE | WHITE | | |
| DECON-CLEAN | WHITE | WHITE | | |
| DECON-QUAT 100 | WHITE | WHITE | | |
| DECON-QUAT 200C | WHITE | WHITE | | |
| DECON-QUAT 200V | WHITE | WHITE | | |
| HYPO-CHLOR 0.25% | WHITE | WHITE | | |
| HYPO-CHLOR 0.52% | WHITE | WHITE | | |
| HYPO-CHLOR 5.25% | WHITE | WHITE | | |
| HYPO-CHLOR Neutral 0.25% | WHITE | WHITE | | |
| HYPO-CHLOR Neutral 0.52% | WHITE | WHITE | | |
| STERI-PEROX 3% | WHITE | WHITE | | |
| STERI-PEROX 6% | WHITE | WHITE | | |
| DECON-SPORE 200 PLUS (SPORICIDE) | WHITE SEMI-TRANSPARENT | WHITE | | |
| DECON-SPORE 200 PLUS (DISINFECTANT) | WHITE SEMI-TRANSPARENT | WHITE | | |
| STEEL-BRIGHT | WHITE | WHITE | | |
| STERI-SILICON | WHITE | BLACK | | |
| DECON-GLASS | WHITE | WHITE | | |
| VAI WFI QUALITY WATER | WHITE | WHITE | | |
| STERI-WATER | WHITE | WHITE | | |

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
 Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
 15JAN2025 / DECON, EPA



Product Labeling

HYPO-CHLOR® Neutral 0.25% and 0.52%

Neutralized Sodium Hypochlorite 0.25% & 0.52% Formulated with Water for Injection

HYPO-CHLOR® NEUTRAL

.25%

Sterile Pharmaceutical Clean Room Formula
Sodium Hypochlorite at 0.25% formulated with USP Water for Injection Quality Water - For Use in Cleanrooms and Controlled Areas - Hard Surface Disinfectant - Fungicidal - Aseptically Filtered at 0.2 µm - Ready to Use After Mixing

Active ingredient:
**Sodium hypochlorite (CAS # 7681-52-8) 0.25%

Other ingredients:
99.75%

Total: 100.00%

**After mixing a decrease in % of active ingredient is expected, use within 7 days or discard any unused material.

KEEP OUT OF REACH OF CHILDREN
WARNING

Net Contents:

- 16 fl. oz. x 12 (Total: 1.5 Gallons)
Reorder # SHC-NPH-0.25-16Z
- 1 Gallon x 4 (Total: 4 Gallons)
Reorder # SHC-NPH-0.25-0Z

EPA Reg. No.: 68959-11
EPA Est. No.: For EPA Est. No. see the final character of the lot #: L = 68959-PA-01, C = 68959-PA-03, T = 68959-PA-02
Label # SHC-NPH-0.25-00 - SDS: VEL-125, VEL-133 - Rev03/24 - # 22
Made in USA
Manufactured By: Veltex Associates, Inc.
15 Lee Blvd., Malvern, PA 19355-1234
Tel: 610-644-8335
www.sterile.com - Patents: www.sterile.com/patents.

FIRST AID
If in Eyes: If splashed in eyes hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. **If on Skin or Clothing:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice. **If Inhaled:** Move person to fresh air. If person is not breathing, call 911 or an ambulance, and then give artificial respiration, preferably mouth-to-mouth if possible. Call a poison control center or doctor for further treatment advice. **If Swallowed:** Call a poison control center or doctor immediately for treatment advice. Have person sip glasses of water if able to swallow. Do not induce vomiting unless told to by a poison control center or doctor. Do not give anything to an unconscious person.

EMERGENCIES: In the case of emergency call: 1-866-928-0788 (toll free) or 1-252-207-0001.

PRECAUTIONARY STATEMENTS
HAZARDS TO HUMANS AND DOMESTIC ANIMALS
WARNING: Causes substantial but temporary eye injury. Causes skin irritation. Harmful if inhaled. Do not breathe spray mist. Do not get into eyes, on skin, or on clothing. Wear protective eyewear (goggles, face shield, or safety glasses), long sleeved shirt and long pants, socks, shoes, and waterproof or chemical-resistant gloves. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse.

ENVIRONMENTAL HAZARDS
Do not apply directly to water, or to areas where surface water is present or to intertidal areas below the mean high water mark. Do not contaminate water when disposing of equipment washwater or rinsate.

PHYSICAL AND CHEMICAL HAZARDS: Contact with acid releases toxic chlorine gas. Do not mix this product with other chemicals.

DIRECTIONS FOR USE
It is a violation of Federal law to use this product in a manner inconsistent with its labeling.
Read the label before using.

This product is not for use on medical devices.

LOCATIONS OF USE: Disinfects clean rooms and controlled areas such as those in hospitals and health care institutions, and biopharmaceutical, pharmaceutical, medical device and diagnostic manufacturing, cosmetic manufacturing, and lab animal facilities. Use on hard non-porous, inanimate surfaces in aseptic filling and gowning rooms, general manufacturing areas and laboratories or on machinery, tools, tables, counters, laminar flow benches, floors, walls, carts, shelves, made of plastic, Plexiglas®, glass, vinyl, glazed porcelain, laminates, glazed tile, and stainless steel. It is compatible with most hard non-porous surface materials. Not for use on food contact surfaces.

PREPARE USE SOLUTION: Follow SimpleMix System Container instructions.

| SIMPLEMIX SYSTEM CONTAINER INSTRUCTIONS: Trigger Spray Bottle label: | SIMPLEMIX SYSTEM CONTAINER INSTRUCTIONS: Gallon Size Bottle label: |
|---|--|
| 1. To prepare use solution, open cap. 2. Peel off inner seal by grasping tab at far edge and pulling off. 3. Firmly push small, inner container all the way down. 4. Replace cap and tighten. 5. Slowly swirl for 15 seconds. 6. Move spray nozzle to open position. 7. Follow directions for use on label. | 1. To prepare use solution, open cap. 2. Peel off inner seal by grasping tab at far edge and pulling off. 3. Firmly push small, inner container all the way down. 4. Replace cap and tighten. 5. Slowly swirl for 15 seconds. 6. Open small side spout and peel off inner seal, as above. 7. Pour solution from small side spout onto surface to be treated or alternate containers. 8. Follow directions for use on label. |

SPORICIDAL: HYPO-CHLOR Neutral .25% is an effective sporicide against *Bacillus subtilis* (ATCC 19665) and *Clostridium sporogenes* (ATCC 3584) on non-porous hard surfaces in 90 minutes. Pre-clean all surfaces. Thoroughly wet surface by immersion or by flooding with the prepared use solution and allow surface to remain visibly wet for a minimum of 90 minutes. Allow surface to air dry or after 30 minutes rise and wipe dry with sterilized cloth wipe if desired. After mixing, use or discard within 7 days.

DISINFECTANT: HYPO-CHLOR Neutral .25% is an effective hard surface disinfectant on inanimate, non-porous, non-food contact surfaces as demonstrated in testing against *Pseudomonas aeruginosa* (ATCC 6442) and *Staphylococcus aureus* (ATCC 6538) in 10 minutes. Pre-clean all surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 10 minutes. Allow surface to air dry or after 10 minutes rise and wipe dry with sterilized cloth wipe if desired. After mixing, use or discard within 7 days.

FUNGICIDAL: HYPO-CHLOR Neutral .25% is an effective fungicide on inanimate, non-porous non-food contact surfaces, demonstrated by the results of ATCC efficacy tests against *Trichophyton interdigitale* (ATCC 5533) in 10 minutes. Pre-clean all surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 10 minutes. Allow surface to air dry, or after 10 minutes, rise and wipe dry with sterilized cloth wipe if desired. After mixing, use or discard within 7 days.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.
Storage: Store in original container in a cool dry area away from direct sunlight and heat to avoid deterioration. In case of spill, flood area with large quantities of water.
Disposal: Product or rinsates that cannot be used must be diluted with water before disposal in a sanitary sewer. Follow Federal/Provincial/State regulations and Local/Municipal ordinances when disposing of this product. Improper disposal of excess product, spray mixture or rinsate is a violation of Federal/Provincial/State Laws. If these wastes cannot be disposed of by use according to label instructions, contact your Federal/Provincial/State or Local/Municipal environmental control agency for guidance.
Container Handling: Non-refillable container. Do not reuse or refill this container to hold materials other than this product. Offer container for recycling if available and dispose of in a sanitary landfill, or by incineration.

Case Label SHC-NPH-0.25-00 Rev03/24

HYPO-CHLOR® NEUTRAL

.52%

Sterile Pharmaceutical Cleanroom Formula
Sodium Hypochlorite at 0.52% formulated with USP Water for Injection Quality Water - For Use in Cleanrooms and Controlled Areas - Hard Surface Disinfectant - Fungicidal - Aseptically Filtered at 0.2 µm - Ready to Use After Mixing

Active ingredient:
**Sodium hypochlorite (CAS # 7681-52-8) 0.52%

Other ingredients:
99.48%

Total: 100.00%

**After mixing a decrease in % of active ingredient is expected, use within 7 days or discard any unused material.

KEEP OUT OF REACH OF CHILDREN
WARNING

Net Contents:

- 16 fl. oz. x 12 (Total: 1.5 Gallons)
Reorder # SHC-NPH-0.52-16Z
- 1 Gallon x 4 (Total: 4 Gallons)
Reorder # SHC-NPH-0.52-0Z

EPA Reg. No.: 68959-9
EPA Est. No.: For EPA Est. No. see the final character of the lot #: L = 68959-PA-01, C = 68959-PA-03, T = 68959-PA-02
Label # SHC-NPH-0.52-00 - SDS: VEL-125, VEL-128 - Rev03/24 - # 22
Made in USA
Manufactured By: Veltex Associates, Inc.
15 Lee Blvd., Malvern, PA 19355-1234
Tel: 610-644-8335
www.sterile.com - Patents: www.sterile.com/patents.

FIRST AID
If in Eyes: If splashed in eyes hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. **If on Skin or Clothing:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice. **If Inhaled:** Move person to fresh air. If person is not breathing, call 911 or an ambulance, and then give artificial respiration, preferably mouth-to-mouth if possible. Call a poison control center or doctor for further treatment advice. **If Swallowed:** Call a poison control center or doctor immediately for treatment advice. Have person sip glasses of water if able to swallow. Do not induce vomiting unless told to by a poison control center or doctor. Do not give anything to an unconscious person.

EMERGENCIES: In the case of emergency call: 1-866-928-0788 (toll free) or 1-252-207-0001.

PRECAUTIONARY STATEMENTS
HAZARDS TO HUMANS AND DOMESTIC ANIMALS
WARNING: Causes substantial but temporary eye injury. Causes skin irritation. Harmful if inhaled. Do not breathe spray mist. Do not get into eyes, on skin, or on clothing. Wear protective eyewear (goggles, face shield, or safety glasses), long sleeved shirt and long pants, socks, shoes, and waterproof or chemical-resistant gloves. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse.

ENVIRONMENTAL HAZARDS
Do not apply directly to water, or to areas where surface water is present or to intertidal areas below the mean high water mark. Do not contaminate water when disposing of equipment washwater or rinsate.

PHYSICAL AND CHEMICAL HAZARDS: Contact with acid releases toxic chlorine gas. Do not mix this product with other chemicals.

DIRECTIONS FOR USE
It is a violation of Federal law to use this product in a manner inconsistent with its labeling.
Read the label before using.

This product is not for use on medical devices.

LOCATIONS OF USE: Disinfects cleanrooms and controlled areas such as those in hospitals and health care institutions, and biopharmaceutical, pharmaceutical, medical device and diagnostic manufacturing, cosmetic manufacturing, and lab animal facilities. Use on hard non-porous, inanimate surfaces in aseptic filling and gowning rooms, general manufacturing areas and laboratories or on machinery, tools, tables, counters, laminar flow benches, floors, walls, carts, shelves, made of plastic, Plexiglas®, glass, vinyl, glazed porcelain, laminates, glazed tile, and stainless steel. It is compatible with most hard non-porous surface materials. Not for use on food contact surfaces.

PREPARE USE SOLUTION: Follow SimpleMix System Container instructions.

| SIMPLEMIX SYSTEM CONTAINER INSTRUCTIONS: Trigger Spray Bottle label: | SIMPLEMIX SYSTEM CONTAINER INSTRUCTIONS: Gallon Size Bottle label: |
|---|--|
| 1. To prepare use solution, open cap. 2. Peel off inner seal by grasping tab at far edge and pulling off. 3. Firmly push small, inner container all the way down. 4. Replace cap and tighten. 5. Slowly swirl for 15 seconds. 6. Move spray nozzle to open position. 7. Follow directions for use on label. | 1. To prepare use solution, open cap. 2. Peel off inner seal by grasping tab at far edge and pulling off. 3. Firmly push small, inner container all the way down. 4. Replace cap and tighten. 5. Slowly swirl for 15 seconds. 6. Open small side spout and peel off inner seal, as above. 7. Pour solution from small side spout onto surface to be treated or alternate containers. 8. Follow directions for use on label. |

SPORICIDAL: HYPO-CHLOR Neutral .52% is an effective sporicide against *Bacillus subtilis* (ATCC 19665) and *Clostridium sporogenes* (ATCC 3584) on non-porous hard surfaces in 90 minutes. Pre-clean all surfaces. Thoroughly wet surface by immersion or by flooding with the prepared use solution and allow surface to remain visibly wet for a minimum of 90 minutes. Allow surface to air dry or after 30 minutes rise and wipe dry with sterilized cloth wipe if desired. After mixing, use or discard within 7 days.

DISINFECTANT: HYPO-CHLOR Neutral .52% is an effective hard surface disinfectant on inanimate, non-porous, non-food contact surfaces as demonstrated in testing against *Pseudomonas aeruginosa* (ATCC 6442) and *Staphylococcus aureus* (ATCC 6538) in 10 minutes. Pre-clean all surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 10 minutes. Allow surface to air dry or after 10 minutes rise and wipe dry with sterilized cloth wipe if desired. After mixing, use or discard within 7 days.

FUNGICIDAL: HYPO-CHLOR Neutral .52% is an effective fungicide on inanimate, non-porous non-food contact surfaces, demonstrated by the results of ATCC efficacy tests against *Trichophyton interdigitale* (ATCC 5533) in 10 minutes. Pre-clean all surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 10 minutes. Allow surface to air dry, or after 10 minutes, rise and wipe dry with sterilized cloth wipe if desired. After mixing, use or discard within 7 days.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.
Storage: Store in original container in a cool dry area away from direct sunlight and heat to avoid deterioration. In case of spill, flood area with large quantities of water.
Disposal: Product or rinsates that cannot be used must be diluted with water before disposal in a sanitary sewer. Follow Federal/Provincial/State regulations and Local/Municipal ordinances when disposing of this product. Improper disposal of excess product, spray mixture or rinsate is a violation of Federal/Provincial/State Laws. If these wastes cannot be disposed of by use according to label instructions, contact your Federal/Provincial/State or Local/Municipal environmental control agency for guidance.
Container Handling: Non-refillable container. Do not reuse or refill this container to hold materials other than this product. Offer container for recycling if available and dispose of in a sanitary landfill, or by incineration.

Case Label SHC-NPH-0.52-00 Rev03/24

(Any specific product label is available upon request.)

Product Labeling (Continued)

HYPO-CHLOR Neutral 0.25%

Sterile Pharmaceutical Cleanroom Formula

Sodium Hypochlorite at 0.25% formulated with USP Water for Injection Quality Water • For Use in Cleanrooms and Controlled Areas • Hard Surface Disinfectant • Fungicidal • Aseptically Filtered at 0.2 µm • Ready to Use After Mixing

Active Ingredient:

**Sodium Hypochlorite.....(CAS # 7681-52-9).....0.25%

Other Ingredients:.....99.75%

Total:.....**100.00%**

**After mixing a decrease in % of active ingredient is expected, use within 7 days or discard any unused material.

KEEP OUT OF THE REACH OF CHILDREN WARNING

EPA Reg. No.: 68959-9

EPA Est. No.: For PA Est. No. see the final character of the lot #: L = 68959-PA-01,

C = 68959-PA-03,T = 68959-PA-02

Label # SHC-NPH-0.25-00

Manufactured by: Veltek Associates, Inc.

15 Lee Blvd. Malvern, PA 19355-1234 USA

Tel: 610-644-8335

www.sterile.com • Patents: www.sterile.com/patents

Made in USA

FIRST AID

If in Eyes: If splashed in eyes hold eye open and rinse slowly and gently with water for 15-20 minutes, Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. **If on Skin or Clothing:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice. **If Inhaled:** Move person to fresh air. If person is not breathing, call 911 or an ambulance, and then give artificial respiration, preferably mouth-to-mouth if possible. Call a poison control center or doctor for further treatment advice. **If Swallowed:** Call a poison control center or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to by a poison control center or doctor. Do not give anything to an unconscious person.

EMERGENCIES: In the case of emergency call: 1-866-928-0789 (toll free) or 1-215-207-0061.

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234

Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com

15JAN2025 / DECON, EPA

Product Labeling (Continued)

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

WARNING. Causes substantial but temporary eye injury. Causes skin irritation. Harmful if inhaled. Do not breathe spray mist. Do not get into eyes, on skin, or on clothing. Wear protective eyewear (goggles, face shield, or safety glasses), long sleeved shirt and long pants, socks, shoes, and waterproof or chemical-resistant gloves. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse.

ENVIRONMENTAL HAZARDS

Do not apply directly to water, or to areas where surface water is present or to intertidal areas below the mean high water mark. Do not contaminate water when disposing of equipment washwater or rinsate.

PHYSICAL AND CHEMICAL HAZARDS: Contact with acid releases toxic chlorine gas. Do not mix this product with other chemicals.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Read the label before using.

This product is not for use on medical devices.

LOCATIONS OF USE: Disinfects clean rooms and controlled areas such as those in hospitals and health care institutions, and biopharmaceutical, pharmaceutical, medical device and diagnostic manufacturing, cosmetic manufacturing, and lab animal facilities. Use on hard non-porous, inanimate, surfaces in aseptic filling and gowning rooms, general manufacturing areas and laboratories or on: machinery, tools, tables, counters, laminar-flow benches, floors, walls, carts, shelves, made of plastic, Plexiglas®, glass, vinyl, glazed porcelain, laminates, glazed tile, and stainless steel. It is compatible with most hard non-porous surface materials. Not for use on food contact surfaces.

PREPARE USE SOLUTION: Follow SimpleMix System Container instructions. For SimpleMix instructions, see page 5-9.

Product Labeling (Continued)

FOR THE FOLLOWING USES, MIX AND USE WITHIN 7 DAYS:

SPORICIDAL: HYPO-CHLOR Neutral .25% is an effective sporicide against *Bacillus subtilis* (ATCC 19659) and *Clostridium sporogenes* (ATCC 3584) on non-porous hard surfaces in 90 minutes. Pre-clean all surfaces. Thoroughly wet surface by immersion or by flooding with the prepared use solution and allow surface to remain visibly wet for a minimum of 90 minutes. Allow surface to air dry, or after 90 minutes rinse and wipe dry with sterilized cloth wiper, if desired. After mixing, use or discard within 7 days.

DISINFECTION: HYPO-CHLOR Neutral .25% is an effective hard surface disinfectant on inanimate, non-porous, non-food contact surfaces as demonstrated in testing against *Pseudomonas aeruginosa* (ATCC 15442), and *Staphylococcus aureus* (ATCC 6538) in 10 minutes. Pre-clean all surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 10 minutes. Allow surface to air dry, or after 10 minutes rinse and wipe dry with sterilized cloth wiper, if desired. After mixing, use or discard within 7 days.

FUNGICIDAL: HYPO-CHLOR Neutral .25% is an effective fungicide on inanimate, non-porous non-food contact surfaces, demonstrated by the results of AOAC Efficacy Tests against *Trichophyton interdigitale* (ATCC 9533) in 10 minutes. Pre-clean all surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 10 minutes. Allow surface to air dry, or after 10 minutes, rinse and wipe dry with sterilized cloth wiper, if desired. After mixing, use or discard within 7 days.

FOR THE FOLLOWING USES, MIX AND USE WITHIN 24 HOURS:

DISINFECTION: HYPO-CHLOR Neutral 0.25% is an effective hard surface disinfectant on inanimate, non-porous non-food contact surfaces demonstrated in testing against *Pseudomonas aeruginosa* (ATCC# 15442), and *Staphylococcus aureus* (ATCC#6538) in 5 minutes. Pre-clean visibly soiled surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 5 minutes. Allow surface to air dry, or after 5 minutes rinse and wipe dry with sterilized cloth wiper, if desired. After mixing, use or discard within 24 hours.

FUNGICIDAL: HYPO-CHLOR Neutral 0.25% is an effective fungicide on inanimate, non-porous non-food contact surfaces, demonstrated by the results of AOAC Efficacy Tests against *Trichophyton interdigitale* (ATCC# 9533) in 5 minutes. Pre-clean visibly soiled surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 5 minutes. Allow surface to air dry, or after 5 minutes, rinse and wipe dry with sterilized cloth wiper, if desired. After mixing, use or discard within 24 hours.

***VIRUCIDAL:** The results of AOAC Efficacy Tests show that HYPO-CHLOR Neutral 0.25% is an effective virucide on inanimate, non-porous, non-food contact surfaces, in the presence of organic load (5% Bovine Serum), demonstrated in testing against *Human Adenovirus 5 (ATCC# VR-5), *Influenza A Virus H3N2 (ATCC# VR-1679), and *Poliovirus Type 1 (ATCC# VR-1562) in 5 minutes. Pre-clean visibly soiled surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 5 minutes. Allow surface to air dry, or after 5 minutes rinse and wipe dry with sterilized cloth wiper, if desired. After mixing, use or discard within 24 hours.

Veltex Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
15JAN2025 / DECON, EPA

Product Labeling (Continued)

FOR THE FOLLOWING USES, MIX AND USE WITHIN 24 HOURS:

SANITIZATION OF NON-FOOD CONTACT SURFACES: *(not for use in California)*

HYPO-CHLOR Neutral 0.25% is an effective sanitizer on inanimate, non-porous non-food contact surfaces demonstrated in testing against *Staphylococcus aureus* (ATCC# 6538) and *Klebsiella pneumoniae* (ATCC# 4352) in 5 minutes. Pre-clean visibly soiled surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 5 minutes.

Allow surface to air dry, or after 5 minutes, rinse and wipe dry with sterilized cloth wiper, if desired. After mixing, use or discard within 24 hours.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.

Storage: Store in original container in a cool dry area away from direct sunlight and heat to avoid deterioration. In case of spill, flood area with large quantities of water.

Disposal: Product or rinsates that cannot be used must be diluted with water before disposal in a sanitary sewer. Follow Federal/Provincial/State regulations and Local/Municipal ordinances when disposing of this product. Improper disposal of excess product, spray mixture or rinsate is a violation of Federal/Provincial/State Laws. If these wastes cannot be disposed of by use according to label instructions, contact your Federal/Provincial/State or Local/Municipal environmental control agency for guidance.

Container Handling: Non-refillable container. Do not reuse or refill this container to hold materials other than this product. Offer container for recycling if available and dispose of in a sanitary landfill, or by incineration.



Product Labeling (Continued)

HYPO-CHLOR Neutral 0.52%

Sterile Pharmaceutical Cleanroom Formula

Sodium Hypochlorite at 0.52% formulated with USP Water for Injection Quality Water • For Use in Cleanrooms and Controlled Areas • Hard Surface Disinfectant • Fungicidal • Aseptically Filtered at 0.2 µm • Ready to Use After Mixing

Active Ingredient:

**Sodium Hypochlorite.....(CAS # 7681-52-9).....0.52%

Other Ingredients:.....99.48%

Total:.....**100.00%**

**After mixing a decrease in % of active ingredient is expected, use within 7 days or discard any unused material.

KEEP OUT OF THE REACH OF CHILDREN WARNING

EPA Reg. No.: 68959-9

EPA Est. No.: For PA Est. No. see the final character of the lot #: L = 68959 –PA-01,
C = 68959-PA-03, T = 68959-PA-02

Label # SHC-NPH-0.52-00

Manufactured by: Veltek Associates, Inc.

15 Lee Blvd. Malvern, PA 19355-1234 USA

Tel: 610-644-8335

www.sterile.com • Patents: www.sterile.com/patents

Made in USA

FIRST AID

If in Eyes: If splashed in eyes hold eye open and rinse slowly and gently with water for 15-20 minutes, Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. **If on Skin or Clothing:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice. **If Inhaled:** Move person to fresh air. If person is not breathing, call 911 or an ambulance, and then give artificial respiration, preferably mouth-to-mouth if possible. Call a poison control center or doctor for further treatment advice. **If Swallowed:** Call a poison control center or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to by a poison control center or doctor. Do not give anything to an unconscious person.

EMERGENCIES: In the case of emergency call: 1-866-928-0789 (toll free) or 1-215-207-0061.

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
15JAN2025 / DECON, EPA

Product Labeling (Continued)

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

WARNING. Causes substantial but temporary eye injury. Causes skin irritation. Harmful if inhaled. Do not breathe spray mist. Do not get into eyes, on skin, or on clothing. Wear protective eyewear (goggles, face shield, or safety glasses), long sleeved shirt and long pants, socks, shoes, and waterproof or chemical-resistant gloves. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse.

ENVIRONMENTAL HAZARDS

Do not apply directly to water, or to areas where surface water is present or to intertidal areas below the mean high water mark. Do not contaminate water when disposing of equipment washwater or rinsate.

PHYSICAL AND CHEMICAL HAZARDS: Contact with acid releases toxic chlorine gas. Do not mix this product with other chemicals.

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Read the label before using.

This product is not for use on medical devices.

LOCATIONS OF USE: Disinfects clean rooms and controlled areas such as those in hospitals and health care institutions, and biopharmaceutical, pharmaceutical, medical device and diagnostic manufacturing, cosmetic manufacturing, and lab animal facilities. Use on hard non-porous, inanimate, surfaces in aseptic filling and gowning rooms, general manufacturing areas and laboratories or on: machinery, tools, tables, counters, laminar-flow benches, floors, walls, carts, shelves, made of plastic, Plexiglas®, glass, vinyl, glazed porcelain, laminates, glazed tile, and stainless steel. It is compatible with most hard non-porous surface materials. Not for use on food contact surfaces.

PREPARE USE SOLUTION: Follow SimpleMix® System Container instructions. For SimpleMix instructions, see page 5-9.

Product Labeling (Continued)

FOR THE FOLLOWING USES, MIX AND USE WITHIN 7 DAYS:

SPORICIDAL: HYPO-CHLOR Neutral .52% is an effective sporicide against *Bacillus subtilis* (ATCC 19659) and *Clostridium sporogenes* (ATCC 3584) on non-porous hard surfaces in 90 minutes. Pre-clean all surfaces. Thoroughly wet surface by immersion or by flooding with the prepared use solution and allow surface to remain visibly wet for a minimum of 90 minutes. Allow surface to air dry, or after 90 minutes rinse and wipe dry with sterilized cloth wiper, if desired. After mixing, use or discard within 7 days.

DISINFECTION: HYPO-CHLOR Neutral .52% is an effective hard surface disinfectant on inanimate, non-porous, non-food contact surfaces as demonstrated in testing against *Pseudomonas aeruginosa* (ATCC 15442), and *Staphylococcus aureus* (ATCC 6538) in 10 minutes. Pre-clean all surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 10 minutes. Allow surface to air dry, or after 10 minutes rinse and wipe dry with sterilized cloth wiper, if desired. After mixing, use or discard within 7 days.

FUNGICIDAL: HYPO-CHLOR Neutral .52% is an effective fungicide on inanimate, non-porous non-food contact surfaces, demonstrated by the results of AOAC Efficacy Tests against *Trichophyton interdigitale* (ATCC 9533) in 10 minutes. Pre-clean all surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 10 minutes. Allow surface to air dry, or after 10 minutes, rinse and wipe dry with sterilized cloth wiper, if desired. After mixing, use or discard within 7 days.

FOR THE FOLLOWING USES, MIX AND USE WITHIN 24 HOURS:

DISINFECTION: HYPO-CHLOR Neutral 0.52% is an effective hard surface disinfectant on inanimate, non-porous non-food contact surfaces demonstrated in testing against *Pseudomonas aeruginosa* (ATCC# 15442), and *Staphylococcus aureus* (ATCC#6538) in 5 minutes. Pre-clean visibly soiled surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 5 minutes. Allow surface to air dry, or after 5 minutes rinse and wipe dry with sterilized cloth wiper, if desired. After mixing, use or discard within 24 hours.

FUNGICIDAL: HYPO-CHLOR Neutral 0.52% is an effective fungicide on inanimate, non-porous non-food contact surfaces, demonstrated by the results of AOAC Efficacy Tests against *Trichophyton interdigitale* (ATCC# 9533) in 5 minutes. Pre-clean visibly soiled surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 5 minutes. Allow surface to air dry, or after 5 minutes, rinse and wipe dry with sterilized cloth wiper, if desired. After mixing, use or discard within 24 hours.

***VIRUCIDAL:** The results of AOAC Efficacy Tests show that HYPO-CHLOR Neutral 0.52% is an effective virucide on inanimate, non-porous, non-food contact surfaces, demonstrated in testing against *Human Adenovirus 5 (ATCC# VR-5), *Influenza A Virus, H3N2 (ATCC# VR-1679), and *Poliovirus Type 1(ATCC# VR-1562) in 5 minutes. Pre-clean visibly soiled surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 5 minutes. Allow surface to air dry, or after 5 minutes rinse and wipe dry with sterilized cloth wiper, if desired. After mixing, use or discard within 24 hours.

Product Labeling (Continued)

FOR THE FOLLOWING USES, MIX AND USE WITHIN 24 HOURS:

SANITIZATION OF NON-FOOD CONTACT SURFACES: *(not for use in California)*

HYPO-CHLOR Neutral 0.52% is an effective sanitizer on inanimate, non-porous non-food contact surfaces demonstrated in testing against *Staphylococcus aureus* (ATCC#6538) and *Klebsiella pneumoniae* (ATCC#4352) in 5 minutes. Pre-clean visibly soiled surfaces. Thoroughly wet surface with the prepared use solution and allow to remain wet for a minimum of 5 minutes.

Allow surface to air dry, or after 5 minutes, rinse and wipe dry with sterilized cloth wiper, if desired. After mixing, use or discard within 24 hours.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.

Storage: Store in original container in a cool dry area away from direct sunlight and heat to avoid deterioration. In case of spill, flood area with large quantities of water.

Disposal: Product or rinsates that cannot be used must be diluted with water before disposal in a sanitary sewer. Follow Federal/Provincial/State regulations and Local/Municipal ordinances when disposing of this product. Improper disposal of excess product, spray mixture or rinsate is a violation of Federal/Provincial/State Laws. If these wastes cannot be disposed of by use according to label instructions, contact your Federal/Provincial/State or Local/Municipal environmental control agency for guidance.

Container Handling: Non-refillable container. Do not reuse or refill this container to hold materials other than this product. Offer container for recycling if available and dispose of in a sanitary landfill, or by incineration.

TESTING SUMMARY

HYPO-CHLOR Neutral Efficacy Summary

HYPO-CHLOR Neutral 0.25% & HYPO-CHLOR Neutral 0.52% have been tested to certified AOAC efficacy standards. The products have achieved passing results supporting bactericidal, fungicidal, and sporicidal claims.

7-day hold: Based on the AOAC certified testing presented below, **HYPO-CHLOR Neutral 0.25% & HYPO-CHLOR Neutral 0.52%** are effective as a bactericide and fungicide at 10 minutes, and a sporicide at 90 minutes.

| AOAC 961.02 | | | | | |
|--|---|---|-------------------|-------------------|-------------------|
| Acceptance Criteria: No visible growth observed in ≤ 3 replicates out of 60 (S. aureus) and in ≤ 6 replicate out of 60 (P. aeruginosa). | | | | | |
| Conditions: Organic Soil Load – N/A Contact Time: 10 minutes Contact Temperature: 20°C ± 1°C Contact Relative Humidity: 48% RH Performed to GLP standards | Test Organism | Carrier Population (CFU/Carrier) | Sample Lot | # Carriers | # Positive |
| | <i>Pseudomonas aeruginosa</i> ATCC #15442 | 3.4 x 10 ⁵ | 1 | 60 | 0/60 |
| | | 2.6 x 10 ⁵ | 2 | 60 | 0/60 |
| | | 2.7 x 10 ⁵ | 3 | 60 | 0/60 |
| | <i>Staphylococcus aureus</i> ATCC # 6538 | 8.8 x 10 ⁵ | 1 | 60 | 0/60 |
| | | 9.0 x 10 ⁵ | 2 | 60 | 0/60 |
| | | 8.4 x 10 ⁵ | 3 | 60 | 0/60 |

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
 Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
 15JAN2025 / DECON, EPA

HYPO-CHLOR Neutral Efficacy Summary (Continued)

| AOAC 961.02 | | | | | |
|---|---|----------------------------------|------------|------------|------------|
| Acceptance Criteria: No visible growth observed in any carrier for each organism. | | | | | |
| Conditions: Organic Soil Load – N/A Contact Time: 10 minutes Contact Temperature: 20°C ± 1°C Contact Relative Humidity: 40-50% RH Performed to GLP standards | Test Organism | Carrier Population (CFU/Carrier) | Sample Lot | # Carriers | # Positive |
| | <i>Trichophyton interdigitale</i> ATCC #9533 | 8.0 x 10 ⁴ | 1 | 10 | 0/10 |
| | | 8.5 x 10 ⁴ | 2 | 10 | 0/10 |

| AOAC 955.15 & 966.04 | | | | | |
|--|---|----------------------------------|------------|------------|------------|
| Acceptance Criteria: No visible growth can be detected in any of the test subculture broth tubes. | | | | | |
| Conditions: Organic Soil Load – N/A Contact Time: 90 minutes Contact Temperature: 20°C ± 1°C Contact Relative Humidity: 26% RH Performed to GLP standards | Test Organism | Carrier Population (CFU/Carrier) | Sample Lot | # Carriers | # Positive |
| | <i>Bacillus subtilis</i> ATCC #19659 | 5.6 x 10 ⁵ | 1 | 60 | 0/60 |
| | | | 2 | 60 | 0/60 |
| | | | 3 | 60 | 0/60 |
| | <i>Clostridium sporogenes</i> ATCC #3548 | 5.1 x 10 ⁵ | 1 | 60 | 0/60 |
| | | | 2 | 60 | 0/60 |
| 3 | | | 60 | 0/60 | |

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
 Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
 15JAN2025 / DECON, EPA

HYPO-CHLOR Neutral Efficacy Summary (Continued)

24-hour hold: Based on the AOAC certified testing presented below, **HYPO-CHLOR Neutral 0.25% & HYPO-CHLOR Neutral 0.52%** are effective as a virucide and bactericide at 5 minutes.

| ASTM E1053 | | | | |
|---|--|-------------------|---|---|
| Acceptance Criteria: Demonstrate at least 3log₁₀ reduction in viral titer. | | | | |
| Conditions: Organic Soil Load – 5.0% Fetal Bovine Serum | Test Organism | Sample Lot | Carrier Population (CFU/Carrier) | Log₁₀ Reduction per Carrier |
| Contact Time: 5 minutes Contact Temperature: 23°C ± 1°C Performed to GLP standards | <i>Human Adenovirus type 5, ATCC VR-5, Strain Adenoid 75</i> | 1 | 7.05 log ₁₀ | ≥6.25 |
| | | 2 | 7.05 log ₁₀ | ≥6.25 |
| | <i>Poliovirus type 1 ATCC VR-1562, Strain Chat</i> | 1 | 6.30 log ₁₀ | ≥5.50 |
| | | 2 | 6.30 log ₁₀ | ≥5.50 |
| | <i>Influenza A (H1N1) ATCC VR-1 469, Strain A/PR/8/34</i> | 1 | 5.80 log ₁₀ | ≥5.00 |
| | | 2 | 5.80 log ₁₀ | ≥5.00 |

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
 Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
 15JAN2025 / DECON, EPA

HYPO-CHLOR Neutral Efficacy Summary (Continued)

| AOAC 961.02 | | | | | | |
|--|--|------------|--------------------|------------------------------------|------------|------|
| Acceptance Criteria: No visible growth observed in ≤ 3 replicates out of 60 (<i>S. aureus</i>) and in ≤ 6 replicate out of 60 (<i>P. aeruginosa</i>). | | | | | | |
| Conditions: Organic Soil Load – 5.0% Fetal Bovine Serum Contact Time: 5 minutes Contact Temperature: 23°C \pm 1°C Performed to GLP standards | Test Organism | Sample Lot | Carrier Population | # Carriers | # Positive | |
| | <i>Pseudomonas aeruginosa</i> ATCC #15442 | A | | 1.46 x 10 ⁵ CFU/Carrier | 60 | 0/60 |
| | | B | | 1.13 x 10 ⁶ CFU/Carrier | 60 | 0/60 |
| | | C | | 1.15 x 10 ⁶ CFU/Carrier | 60 | 0/60 |
| | <i>Staphylococcus aureus</i> ATCC # 6538 | A | | 2.40 x 10 ⁶ CFU/Carrier | 60 | 0/60 |
| | | B | | 1.26 x 10 ⁶ CFU/Carrier | 60 | 0/60 |
| | | C | | 1.33 x 10 ⁶ CFU/Carrier | 60 | 0/60 |

| AOAC 961.02 | | | | | | |
|--|---|------------|--------------------|------------------------------------|------------|------|
| Acceptance Criteria: No visible growth observed in any carrier for each organism. | | | | | | |
| Conditions: Organic Soil Load – 5.0% Fetal Bovine Serum Contact Time: 5 minutes Contact Temperature: 24°C \pm 1°C Performed to GLP standards | Test Organism | Sample Lot | Carrier Population | # Carriers | # Positive | |
| | <i>Trichophyton interdigitale</i> ATCC #9533 | A | | 4.00 x 10 ⁴ CFU/Carrier | 10 | 0/10 |
| | | B | | 4.00 x 10 ⁴ CFU/Carrier | 10 | 0/10 |

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
 Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
 15JAN2025 / DECON, EPA

HYPO-CHLOR Neutral Efficacy Summary (Continued)

| ASTN E1153 | | | | | |
|--|---|-----------------------------------|-------------------|-------------------|-------------------|
| Acceptance Criteria: The test substance must demonstrate a minimum 99.9% reduction of test survivors as compared to the population control. | | | | | |
| Conditions: Organic Soil Load – N/A | Test Organism | Carrier Population | Sample Lot | # Carriers | # Positive |
| Contact Time: 5 minutes Contact Temperature: 20°C ± 1°C Performed to GLP standards | <i>Enterobacter aerogenes</i> ATCC #13048 | 2.57×10^6 CFU/Carrier | A | 5 | >99.9% |
| | | | B | 5 | >99.9% |
| | | | C | 5 | >99.9% |
| | <i>Staphylococcus aureus</i> ATCC #6538 | 4.37×10^6 CFU/Carrier | A | 5 | >99.9% |
| | | | B | 5 | >99.9% |
| | | | C | 5 | >99.9% |

ADDITIONAL DOCUMENTATION

Upon request, the following additional documentation is available:

- Product Testing Reports
- Product Labels
- Product Validation
- Safety Data Sheet
- Sample Lot Specific Documentation - Delivered with lot specific analytical and sterility data, tested to current USP compendium



VAI®'s SCMD Product Family

VAI's Sterile Chemical Manufacturing Division (SCMD) manufactures a complete range of cleaning agents and disinfectants that are used daily in cleanroom operations. Overall, VAI's capabilities for manufacturing products include the ability to fill aerosol, bulk, and unit dose packages in ISO 5 or 7 (Grade A/B). Our aseptic filling operations are coupled with the validated and proven ability to irradiate a final product. Assurances are taken in every aspect of SCMD concerning sterility and particulate removal. VAI's operations mirror current GMP's and enforces the adherence to USP specifications. VAI is an EPA and FDA registered facility.

For more information call 610-644-8335, email sales@sterile.com or visit our website at www.sterile.com. Patents: www.sterile.com/patents

Veltek Associates, Inc.

15 Lee Boulevard, Malvern, PA 19355-1234
Tel: 610-644-8335 Fax: 610-644-8335 www.sterile.com
15JAN2025 / DECON, EPA