

**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1. Product identifier**

Product form : Mixture
Product name : DEC-CYCLE®. Also known as: DEC-CYCLE® II
Product code : SDS DCY-98-01-EU

1.2. Relevant identified uses of the substance or mixture and uses advised against**1.2.1. Relevant identified uses**

Use of the substance/mixture : Disinfectant
Product for industrial use only

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Veltek Associates, Inc.
15 Lee Blvd
Malvern, PA 19355-1234 USA
Telephone: +1 610-644-8335 - Fax: +1 610-644-8336
E-mail: vai@sterile.com

Veltek Associates Inc., Branch Office Europe
PO Box 1062, 8200 BB Lelystad, Netherlands
Customer service (USA): +800 00888700

India distributor:
M/s. Shah Brothers
C-32, Shri Ram Indl. Estate
G.D. Ambekar Marg
Wadala, Mumbai- 400031 India
Telephone: +91 22-43560400

1.4. Emergency telephone number

Emergency number : For Spill/Exposure Emergency Response Service in Europe in English (and 23 other European languages) (24 hours): +44 (0)1235 239 670
Ireland Poison Centre: +353 (0)1 809 2166 (Available to the public 08.00–22.00) (Ireland only)
For Middle East/Africa (24 hours): +44 (0)1235 239 671
For Hindi (24 hours): 000 800 100 7479

SECTION 2: Hazards identification**2.1. Classification of the substance or mixture****Classification according to Regulation (EC) No. 1272/2008 [CLP]**

Flam. Liq. 3	H226
Skin Irrit. 2	H315
Eye Dam. 1	H318
Skin Sens. 1	H317
Carc. 2	H351
Repr. 2	H361f
Aquatic Chronic 1	H410

Full text of hazard classes and H-statements : see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

DEC-CYCLE®. Also known as: DEC-CYCLE® II

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)



Signal word (CLP)

Hazardous ingredients

Hazard statements (CLP)

Precautionary statements (CLP)

- : Danger
- : Biphenyl-2-ol; Clorofene
- : H226 - Flammable liquid and vapour.
H315 - Causes skin irritation.
H317 - May cause an allergic skin reaction.
H318 - Causes serious eye damage.
H351 - Suspected of causing cancer.
H361f - Suspected of damaging fertility.
H410 - Very toxic to aquatic life with long lasting effects.
- : P201 - Obtain special instructions before use.
P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P280 - Wear eye protection, protective gloves, protective clothing.
P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310 - Immediately call a doctor.
P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

2.3. Other hazards

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII

This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
2-methylpentane-2,4-diol	(CAS-No.) 107-41-5 (EC No.) 203-489-0 (EC index No.) 603-053-00-3	10 - < 30	Eye Irrit. 2, H319 Skin Irrit. 2, H315
Citric acid	(CAS-No.) 77-92-9 (EC No.) 201-069-1	10 - < 30	Eye Irrit. 2, H319
Isopropanol	(CAS-No.) 67-63-0 (EC No.) 200-661-7 (EC index No.) 603-117-00-0	7 - < 13	Flam. Liq. 2, H225 Eye Irrit. 2, H319 STOT SE 3, H336
Biphenyl-2-ol	(CAS-No.) 90-43-7 (EC No.) 201-993-5 (EC index No.) 604-020-00-6	5 - < 10	Skin Irrit. 2, H315 Eye Dam. 1, H318 STOT SE 3, H335 Aquatic Acute 1, H400 Aquatic Chronic 1, H410

DEC-CYCLE®. Also known as: DEC-CYCLE® II

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

Clorofene	(CAS-No.) 120-32-1 (EC No.) 204-385-8 (EC index No.) 604-093-00-4	5 - < 10	Acute Tox. 4 (Inhalation), H332 Skin Irrit. 2, H315 Eye Dam. 1, H318 Skin Sens. 1, H317 Carc. 2, H351 Repr. 2, H361f STOT RE 2, H373 Aquatic Acute 1, H400 Aquatic Chronic 1, H410 (M=100)
2,2',2"-nitrotriethanol	(CAS-No.) 102-71-6 (EC No.) 203-049-8	1 - < 5	Not classified
Benzenesulfonic acid, C10-16-alkyl derivs.	(CAS-No.) 68584-22-5 (EC No.) 271-528-9	1 - < 5	Eye Irrit. 2, H319

Full text of H-statements: see section 16

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general	: Never give anything by mouth to an unconscious person. If you feel unwell, seek medical advice (show the label where possible).
First-aid measures after inhalation	: Remove person to fresh air and keep comfortable for breathing. If symptoms develop obtain medical attention.
First-aid measures after skin contact	: Remove affected clothing and wash all exposed skin area with mild soap and water, followed by warm water rinse. If skin irritation or rash occurs: Get medical advice/attention.
First-aid measures after eye contact	: Rinse immediately with plenty of water (for at least 15 minutes). Ensure that folded skin of eyelids is thoroughly washed with water. Remove contact lenses, if present and easy to do. Continue rinsing. Obtain immediate medical attention.
First-aid measures after ingestion	: Do NOT induce vomiting. Do not give an unconscious person anything to drink. Wash out mouth with water and give 100-200 ml of water to drink. If symptoms develop, obtain medical attention.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects after inhalation	: At high concentrations, the vapours can be irritating to the respiratory system.
Symptoms/effects after skin contact	: Causes skin irritation. May cause an allergic skin reaction. Skin rash/inflammation.
Symptoms/effects after eye contact	: Causes serious eye damage.
Symptoms/effects after ingestion	: Ingestion may cause irritation of the gastrointestinal tract.
Chronic symptoms	: Suspected of causing cancer. Suspected of damaging fertility or the unborn child.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media	: Alcohol-resistant foam. Dry powder. Carbon dioxide. Water spray.
Unsuitable extinguishing media	: Do not use a heavy water stream.

5.2. Special hazards arising from the substance or mixture

Fire hazard	: Flammable liquid and vapour. Vapours are heavier than air and may travel considerable distance to an ignition source and flash back to source of vapours.
Explosion hazard	: Containers may rupture when heated.
Hazardous decomposition products in case of fire	: Carbon dioxide. Carbon monoxide. Chlorine. Hydrocarbons.

5.3. Advice for firefighters

Firefighting instructions	: Exercise caution when fighting any chemical fire. Use water spray or fog for cooling exposed containers. Avoid fire-fighting water entering the environment.
---------------------------	--

DEC-CYCLE®. Also known as: DEC-CYCLE® II

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection. Use self-contained breathing apparatus when in close proximity to fire.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Emergency procedures : Remove all sources of ignition. Ventilate area. Do not breathe vapours. Do not get in eyes, on skin, or on clothing. Evacuate unnecessary personnel.

6.1.2. For emergency responders

Protective equipment : Equip cleanup crew with proper protection. Use chemically protective clothing.
Emergency procedures : Remove ignition sources. Use only non-sparking tools. Take precautionary measures against static discharge. Ventilate area. Do not breathe vapours. Do not get in eyes, on skin, or on clothing.

6.2. Environmental precautions

Collect spillage. Prevent entry to sewers and public waters. Notify authorities if large amounts of the product enters sewers or public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Use non-sparking tools. Absorb with earth, sand or other non-combustible material and transfer to containers for later disposal.

6.4. Reference to other sections

SECTION 8: Exposure controls/personal protection. SECTION 13: Disposal considerations.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Ground/bond container and receiving equipment. Take precautionary measures against static discharge. Use only non-sparking tools. Use only outdoors or in a well-ventilated area. Provide good ventilation in process area to prevent formation of vapour. Do not breathe vapours. Do not get in eyes, on skin, or on clothing.
Hygiene measures : Handle in accordance with good industrial hygiene and safety practice. Do not eat, drink or smoke when using this product. Take off immediately all contaminated clothing and wash it before reuse. Contaminated work clothing should not be allowed out of the workplace. Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Protect material from direct sunlight. Store in original container. Keep container tightly closed. Store in a well-ventilated place. Keep cool. Store locked up.
Incompatible materials : Oxidising agents. Reducing agents. Strong alkalis. Strong acids.

7.3. Specific end use(s)

Product for industrial use only. Disinfectant.

DEC-CYCLE®. Also known as: DEC-CYCLE® II

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

2-methylpentane-2,4-diol (107-41-5)

Ireland - Occupational Exposure Limits

Local name	Hexylene glycol
OEL (15 min ref) (mg/m ³)	125 mg/m ³
OEL (15 min ref) (ppm)	25 ppm
Regulatory reference	Code of Practice for the Chemical Agents Regulations 2018

United Kingdom - Occupational Exposure Limits

Local name	2-Methylpentane-2,4-diol
WEL TWA (mg/m ³)	123 mg/m ³
WEL TWA (ppm)	25 ppm
WEL STEL (mg/m ³)	123 mg/m ³
WEL STEL (ppm)	25 ppm
Regulatory reference	EH40/2005 (Third edition, 2018). HSE

Isopropanol (67-63-0)

Ireland - Occupational Exposure Limits

Local name	Isopropyl alcohol
OEL (8 hours ref) (ppm)	200 ppm
OEL (15 min ref) (ppm)	400 ppm
Notes (IE)	Sk
Regulatory reference	Code of Practice for the Chemical Agents Regulations 2018

United Kingdom - Occupational Exposure Limits

Local name	Propan-2-ol
WEL TWA (mg/m ³)	999 mg/m ³
WEL TWA (ppm)	400 ppm
WEL STEL (mg/m ³)	1250 mg/m ³
WEL STEL (ppm)	500 ppm
Regulatory reference	EH40/2005 (Third edition, 2018). HSE

2,2',2''-nitrilotriethanol (102-71-6)

Ireland - Occupational Exposure Limits

Local name	Triethanolamine
OEL (8 hours ref) (mg/m ³)	5 mg/m ³
Regulatory reference	Code of Practice for the Chemical Agents Regulations 2018

8.2. Exposure controls

Appropriate engineering controls:

Provide good ventilation in process area to prevent formation of vapour. Ensure exposure is below occupational exposure limits (where available). Local exhaust ventilation (LEV) may be required to control inhalation exposure. Emergency eye wash fountains should be available in the immediate vicinity of any potential exposure.

DEC-CYCLE®. Also known as: DEC-CYCLE® II

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

Personal protective equipment:

Avoid all unnecessary exposure.

Hand protection:

Wear suitable gloves resistant to chemical penetration. Standard EN 374 - Protective gloves against chemicals. The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed. Gloves should be removed and replaced if there are any signs of degradation or breakthrough.

Eye protection:

Chemical goggles or safety glasses. Standard EN 166 - Personal eye-protection.

Skin and body protection:

Wear chemically protective gloves, lab coat or apron to prevent prolonged or repeated skin contact

Respiratory protection:

In case of inadequate ventilation wear respiratory protection. EN 405

Thermal hazard protection:

Not required for normal conditions of use.

Environmental exposure controls:

Avoid release to the environment.

Other information:

Do not eat, drink or smoke during use.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	: Liquid
Appearance	: Clear.
Colour	: Amber.
Odour	: Alcohol.
Odour threshold	: No data available
pH	: 1.5 – 2.5 (Concentrate)
Relative evaporation rate (butylacetate=1)	: No data available
Melting point	: No data available
Freezing point	: No data available
Boiling point	: 100 °C
Flash point	: 41.6 °C
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Flammability (solid, gas)	: Not applicable
Vapour pressure	: No data available
Relative vapour density at 20 °C	: No data available
Relative density	: 1.092 +/-0.01 (25 °C) (Water = 1)
Solubility	: Water: Miscible
Log Pow	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: Vapours may form explosive mixture with air.
Oxidising properties	: Not oxidising.
Explosive limits	: No data available

9.2. Other information

No additional information available

DEC-CYCLE®. Also known as: DEC-CYCLE® II

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under recommended handling and storage conditions (see section 7).

10.2. Chemical stability

Stable under recommended handling and storage conditions (see section 7).

10.3. Possibility of hazardous reactions

Under fire conditions closed containers may rupture or explode.

10.4. Conditions to avoid

Keep away from open flames, hot surfaces and sources of ignition. Extremely high or low temperatures. Direct sunlight.

10.5. Incompatible materials

Oxidising agents. Reducing agents. Strong alkalis. Strong acids.

10.6. Hazardous decomposition products

Fire may produce irritating, corrosive and/or toxic gases. Carbon monoxide. Carbon dioxide. Chlorine. Hydrocarbons.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity (oral)	: Not classified
Acute toxicity (dermal)	: Not classified
Acute toxicity (inhalation)	: Not classified
Additional information	: Based on available data, the classification criteria are not met

2-methylpentane-2,4-diol (107-41-5)

LD50 oral, rat	> 2000 mg/kg
LD50 dermal, rabbit	> 2000 mg/kg
LC50 inhalation, rat (mg/l)	> 60 ml/m ³ - 8 Hours

Isopropanol (67-63-0)

LD50 oral, rat	5840 mg/kg
LD50 dermal, rat	16.4 ml/kg
LC50 inhalation, rat (ppm)	> 10000 ppm - 6 Hours

Citric acid (77-92-9)

LD50 oral, rat	5400 mg/kg
LD50 dermal, rat	> 2000 mg/kg

Biphenyl-2-ol (90-43-7)

LD50 oral, rat	2733 mg/kg
LD50 dermal, rabbit	> 2000 mg/kg
LC50 inhalation, rat (mg/l)	> 949 mg/m ³ - 1 Hours

Clorofene (120-32-1)

LD50 oral, rat	4147 mg/kg (OECD 401 method)
----------------	------------------------------

DEC-CYCLE®. Also known as: DEC-CYCLE® II

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

LD50 dermal, rat	> 2000 mg/kg (OECD 402 method)
LC50 inhalation, rat (mg/l)	2.5 mg/l - 4 Hours (OECD 403 method)

2,2',2"-nitrioltriethanol (102-71-6)

LD50 oral, rat	6400 mg/kg (OECD 401 method)
LD50 dermal, rat	> 2000 mg/kg (OECD 402 method)

Skin corrosion/irritation	: Causes skin irritation. pH: 1.5 – 2.5 (Concentrate)
Serious eye damage/irritation	: Causes serious eye damage. pH: 1.5 – 2.5 (Concentrate)
Respiratory or skin sensitisation	: May cause an allergic skin reaction.
Germ cell mutagenicity	: Not classified
Additional information	: Based on available data, the classification criteria are not met
Carcinogenicity	: Suspected of causing cancer.

Isopropanol (67-63-0)

IARC group	3 - Not classifiable
------------	----------------------

Biphenyl-2-ol (90-43-7)

IARC group	3 - Not classifiable
------------	----------------------

2,2',2"-nitrioltriethanol (102-71-6)

IARC group	3 - Not classifiable
------------	----------------------

Reproductive toxicity	: Suspected of damaging fertility.
STOT-single exposure	: Not classified
Additional information	: Based on available data, the classification criteria are not met
STOT-repeated exposure	: Not classified
Additional information	: Based on available data, the classification criteria are not met
Aspiration hazard	: Not classified
Additional information	: Based on available data, the classification criteria are not met
Potential adverse human health effects and symptoms	: Causes serious eye damage. Causes skin irritation. May cause an allergic skin reaction. Skin rash/inflammation. At high concentrations, the vapours can be irritating to the respiratory system. Ingestion may cause irritation of the gastrointestinal tract. Suspected of damaging fertility. Suspected of causing cancer.

SECTION 12: Ecological information

12.1. Toxicity

Hazardous to the aquatic environment, short-term (acute)	: Not classified
Hazardous to the aquatic environment, long-term (chronic)	: Very toxic to aquatic life with long lasting effects.

2-methylpentane-2,4-diol (107-41-5)

LC50 fish	10000 mg/l 96 Hours - Lepomis macrochirus
EC50 Daphnia	2700 – 3700 mg/l 48 Hours - Daphnia magna

DEC-CYCLE®. Also known as: DEC-CYCLE® II

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

Isopropanol (67-63-0)	
LC50 fish	9640 mg/l - 96 Hours (Pimephales promelas)
EC50 Daphnia	> 10000 mg/l - 48 Hours (Daphnia magna)
NOEC chronic algae	1800 mg/l - 7 days (Scenedesmus quadricauda)

Citric acid (77-92-9)	
LC50 fish	440 – 760 mg/l 48 Hours - Leuciscus idus melanotus
EC50 Daphnia	1535 mg/l 24 Hours - Daphnia magna

Biphenyl-2-ol (90-43-7)	
LC50 fish	4.5 mg/l - 96 Hours (Danio rerio)
EC50 Daphnia	2.7 mg/l - 48 Hours (Daphnia magna, Mobility)
ErC50 (algae)	3.57 mg/l - 72 Hours (Pseudokirchneriella subcapitata, Growth rate), (OECD 201 method)
NOEC (chronic)	0.009 mg/l - 21 days (Daphnia magna, reproduction), (OECD 211 method)
NOEC chronic fish	0.036 mg/l - 21 days (Pimephales promelas, reproduction)
NOEC chronic algae	0.468 mg/l - 72 Hours (Pseudokirchneriella subcapitata, Growth rate), (OECD 201 method)

Clorofene (120-32-1)	
LC50 fish	1.5 mg/l - 96 Hours (Danio rerio)
EC50 Daphnia	0.655 mg/l - 48 Hours (Daphnia magna)(OECD 202 method)
EC50 other aquatic organisms 1	0.089 mg/l - 96 Hours (Americamysis bahia)(EPA OPPTS 850.1035)
ErC50 (algae)	0.435 mg/l - 72 Hours (Navicula pelliculosa, Growth rate)(OCSP 850.4500)
ErC50 (other aquatic plants)	0.155 mg/l - 96 Hours (Skeletonema costatum, Growth rate)((OECD 201 method)
NOEC chronic fish	< 0.0095 mg/l - 30 days (Danio rerio) (OECD 210 method)
NOEC chronic crustacea	0.0067 mg/l - 21 days (Daphnia magna, reproduction)

2,2',2''-nitrioltriethanol (102-71-6)	
LC50 fish	11800 mg/l - 96 Hours (Pimephales promelas)(APHA)
EC50 Daphnia	609.88 mg/l - 96 Hours (Ceriodaphnia dubia)(ASTM E1192)
EC50 72h algae (1)	216 mg/l - 72 Hours (Desmodesmus subspicatus, Growth rate)(DIN 38412, 9)
NOEC (acute)	16 mg/l 21 days - Daphnia magna
NOEC chronic crustacea	16 mg/l - 21 days (Daphnia magna)

12.2. Persistence and degradability

DEC-CYCLE®. Also known as: DEC-CYCLE® II	
Persistence and degradability	Readily biodegradable.
Biodegradation	74 % - 28 days

Isopropanol (67-63-0)	
Persistence and degradability	Expected to be readily biodegradable.

DEC-CYCLE®. Also known as: DEC-CYCLE® II

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

Biphenyl-2-ol (90-43-7)	
Persistence and degradability	Readily biodegradable.
Biodegradation	70.8 – 75.7 % - 28 days (OECD 301B method)

Clorofene (120-32-1)	
Persistence and degradability	Inherently biodegradable. Not readily biodegradable.

12.3. Bioaccumulative potential

DEC-CYCLE®. Also known as: DEC-CYCLE® II	
Bioaccumulative potential	No information available.

Isopropanol (67-63-0)	
Bioconcentration factor (BCF REACH)	3
Log Pow	0.05

Biphenyl-2-ol (90-43-7)	
BCF fish 1	21.7 (Danio rerio)
Log Pow	3.18 (22.5 °C)(OECD 107 method)
Bioaccumulative potential	Not expected to bioaccumulate.

Clorofene (120-32-1)	
Log Pow	4.276 (25 °C)
Bioaccumulative potential	Not expected to bioaccumulate.

12.4. Mobility in soil

DEC-CYCLE®. Also known as: DEC-CYCLE® II	
Ecology - soil	Miscible with water.

Biphenyl-2-ol (90-43-7)	
Log Koc	2.4 – 2.6 (20 °C)

12.5. Results of PBT and vPvB assessment

DEC-CYCLE®. Also known as: DEC-CYCLE® II	
This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII	
This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII	

12.6. Other adverse effects

Additional information : Avoid release to the environment.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste disposal recommendations : Do not discharge into drains or the environment. Dispose in a safe manner in accordance with local/national regulations.

Additional information : Handle empty containers with care because residual vapours are flammable. Empty containers should be taken for recycling, recovery or waste in accordance with local regulation.

DEC-CYCLE®. Also known as: DEC-CYCLE® II

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

Ecology - waste materials

: Avoid release to the environment.

SECTION 14: Transport information

In accordance with ADR / IMDG / IATA

14.1. UN number

UN-No. (ADR) : UN 2924
UN-No. (IMDG) : UN 2924
UN-No. (IATA) : UN 2924


14.2. UN proper shipping name

Proper Shipping Name : FLAMMABLE LIQUID, CORROSIVE, N.O.S. (Isopropanol ; Citric acid)
Proper Shipping Name (IMDG) : FLAMMABLE LIQUID, CORROSIVE, N.O.S. (Isopropanol ; Citric acid)
Proper Shipping Name (IATA) : Flammable liquid, corrosive, n.o.s. (Isopropanol ; Citric acid)
Transport document description (ADR) : UN 2924 FLAMMABLE LIQUID, CORROSIVE, N.O.S. (Isopropanol ; Citric acid), 3 (8), III, (D/E), ENVIRONMENTALLY HAZARDOUS
Transport document description (IMDG) : UN 2924 FLAMMABLE LIQUID, CORROSIVE, N.O.S. (Isopropanol ; Citric acid), 3 (8), III, MARINE POLLUTANT/ENVIRONMENTALLY HAZARDOUS
Transport document description (IATA) : UN 2924 Flammable liquid, corrosive, n.o.s. (Isopropanol ; Citric acid), 3 (8), III, ENVIRONMENTALLY HAZARDOUS

14.3. Transport hazard class(es)

ADR

Transport hazard class(es) (ADR) : 3 (8)
Hazard labels : 3, 8




IMDG

Transport hazard class(es) (IMDG) : 3 (8)
Danger labels (IMDG) : 3, 8



IATA

Transport hazard class(es) (IATA) : 3 (8)
Danger labels (IATA) : 3, 8



14.4. Packing group

Packing group : III
Packing group (IMDG) : III
Packing group (IATA) : III

14.5. Environmental hazards

Dangerous for the environment : Yes
Marine pollutant : Yes

DEC-CYCLE®. Also known as: DEC-CYCLE® II

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

Other information : No supplementary information available

14.6. Special precautions for user

Special transport precautions : ADR: 5.2.1.8.1: Environmentally hazardous substances derogation applies (quantity of liquids ≤ 5 litres or net mass of solids ≤ 5 kg)

IMDG 2.10.2.7: Environmentally hazardous substances derogation applies (quantity of liquids ≤ 5 litres or net mass of solids ≤ 5 kg)

Overland transport

Tunnel restriction code (ADR) : D/E

Transport by sea

No data available

Air transport

No data available

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Authorisations and/or restrictions on use (Annex XVII):		
Reference code	Applicable on	Entry title or description
3.	DEC-CYCLE®. Also known as: DEC-CYCLE® II ; 2-methylpentane-2,4-diol ; Isopropanol ; Benzenesulfonic acid, C10-16-alkyl derivs.	Liquid substances or mixtures which are regarded as dangerous in accordance with Directive 1999/45/EC or are fulfilling the criteria for any of the following hazard classes or categories set out in Annex I to Regulation (EC) No 1272/2008
3(a)	DEC-CYCLE®. Also known as: DEC-CYCLE® II ; Isopropanol	Substances or mixtures fulfilling the criteria for any of the following hazard classes or categories set out in Annex I to Regulation (EC) No 1272/2008: Hazard classes 2.1 to 2.4, 2.6 and 2.7, 2.8 types A and B, 2.9, 2.10, 2.12, 2.13 categories 1 and 2, 2.14 categories 1 and 2, 2.15 types A to F
3(b)	DEC-CYCLE®. Also known as: DEC-CYCLE® II ; 2-methylpentane-2,4-diol ; Isopropanol ; Benzenesulfonic acid, C10-16-alkyl derivs.	Substances or mixtures fulfilling the criteria for any of the following hazard classes or categories set out in Annex I to Regulation (EC) No 1272/2008: Hazard classes 3.1 to 3.6, 3.7 adverse effects on sexual function and fertility or on development, 3.8 effects other than narcotic effects, 3.9 and 3.10
3(c)	DEC-CYCLE®. Also known as: DEC-CYCLE® II	Substances or mixtures fulfilling the criteria for any of the following hazard classes or categories set out in Annex I to Regulation (EC) No 1272/2008: Hazard class 4.1
40.	DEC-CYCLE®. Also known as: DEC-CYCLE® II ; Isopropanol	Substances classified as flammable gases category 1 or 2, flammable liquids categories 1, 2 or 3, flammable solids category 1 or 2, substances and mixtures which, in contact with water, emit flammable gases, category 1, 2 or 3, pyrophoric liquids category 1 or pyrophoric solids category 1, regardless of whether they appear in Part 3 of Annex VI to Regulation (EC) No 1272/2008 or not.

Contains no substance on the REACH candidate list

Contains no REACH Annex XIV substances

Contains no substance subject to Regulation (EU) No 649/2012 of the European Parliament and of the Council of 4 July 2012 concerning the export and import of hazardous chemicals.

Contains no substance subject to Regulation (EU) No 2019/1021 of the European Parliament and of the Council of 20 June 2019 on persistent organic pollutants

15.1.2. National regulations

No additional information available

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

DEC-CYCLE®. Also known as: DEC-CYCLE® II

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

SECTION 16: Other information

Indication of changes:			
Section	Changed item	Change	Comments
2	Hazards identification	Modified	
3	Composition/information on ingredients	Modified	
4	First aid measures	Modified	
6	Accidental release measures	Modified	
7	Handling and storage	Modified	
11	Toxicological information	Modified	
12.	Ecological information	Modified	
14	Transport information	Modified	
15	Regulatory information	Modified	

Abbreviations and acronyms:	
	ADR (Accord européen relatif au transport international des marchandises Dangereuses par Route)
	ATE (Acute Toxicity Estimate)
	CAS (Chemical Abstracts Service) number
	CLP (Classification, Labeling and Packaging)
	DNEL (Derived No Effect Level)
	EC (European Community)
	EC50 (Effective Concentration 50%)
	EN (European Norm)
	IARC (International Agency for Research on Cancer)
	IATA (International Air Transport Association)
	IMDG (International Maritime Dangerous Goods Code)
	IMO (International Maritime Organisation)
	LC50 (Lethal Concentration 50%)
	LD50 (Lethal Dose 50%)
	MAC (Maximal Allowed Concentration)
	OECD (Organisation for Economic Co-operation and Development)
	PBT (Persistent, Bioaccumulative and Toxic)
	PNEC (Predicted No Effect Concentration)
	REACH (Registration, Evaluation and Authorisation of CHemicals)
	RID (Règlement concernant le transport international ferroviaire de marchandises)
	STEL (Short Term Exposure Limit)
	TWA (Time Weighted Average)
	UNxxxx (Number assigned by the United Nations Committee of Experts on the Transport of Dangerous Goods)
	vPvB (very Persistent and very Bioaccumulative)

DEC-CYCLE®. Also known as: DEC-CYCLE® II

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

Data sources	: REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.
Other information	: Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]. Physical hazards: On the basis of test data. Health hazards: Calculation method. Environmental hazards: Calculation method.

Classification according to Regulation (EC) No. 1272/2008 [CLP]:	
Flam. Liq. 3	H226
Skin Irrit. 2	H315
Eye Dam. 1	H318
Skin Sens. 1	H317
Carc. 2	H351
Repr. 2	H361f
Aquatic Chronic 1	H410

Full text of H- and EUH-statements:	
Acute Tox. 4 (Inhalation)	Acute toxicity (inhal.), Category 4
Aquatic Acute 1	Hazardous to the aquatic environment — Acute Hazard, Category 1
Aquatic Chronic 1	Hazardous to the aquatic environment — Chronic Hazard, Category 1
Carc. 2	Carcinogenicity, Category 2
Eye Dam. 1	Serious eye damage/eye irritation, Category 1
Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
Flam. Liq. 2	Flammable liquids, Category 2
Flam. Liq. 3	Flammable liquids, Category 3
Repr. 2	Reproductive toxicity, Category 2
Skin Irrit. 2	Skin corrosion/irritation, Category 2
Skin Sens. 1	Skin sensitisation, Category 1
STOT RE 2	Specific target organ toxicity — Repeated exposure, Category 2
STOT SE 3	Specific target organ toxicity — Single exposure, Category 3, Narcosis
STOT SE 3	Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation
H225	Highly flammable liquid and vapour.
H226	Flammable liquid and vapour.
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H318	Causes serious eye damage.
H319	Causes serious eye irritation.
H332	Harmful if inhaled.
H335	May cause respiratory irritation.
H336	May cause drowsiness or dizziness.
H351	Suspected of causing cancer.
H361f	Suspected of damaging fertility.
H373	May cause damage to organs through prolonged or repeated exposure.

DEC-CYCLE®. Also known as: DEC-CYCLE® II

Safety Data Sheet

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2015/830

H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects.

SDS EU (REACH Annex II)

Information contained in this publication or as otherwise supplied to Users is believed to be accurate and is given in good faith, but it is for the Users to satisfy themselves of the suitability of the product for their own particular purpose. Veltek Associates, Inc. gives no warranty as to the fitness of the product for any particular purpose and any implied warranty or condition (statutory or otherwise) is excluded except to the extent that exclusion is prevented by law. Veltek Associates, Inc. accepts no liability for loss or damage resulting from reliance on this information. Freedom under Patents, Copyright and Designs cannot be assumed.

This SDS has been translated into the official language of the country/region in which the product is to be placed on the market. Where no official translation exists, the regulatory text is reported in English, as it appears in the relevant regulatory text.