



DECON-CLEAN[®] WFI FORMULA AEROSOL

Safety Data Sheet

VELTEK ASSOCIATES, INC.

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations
Date of issue: 10/05/2011 Revision date: 01/26/2017 Supersedes: 04/20/2016 Version: 3.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Mixture
Product name : DECON-CLEAN[®] WFI FORMULA AEROSOL
Product code : SDS VEL-031

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/mixture : Residue Remover
Industrial/Professional use spec : For professional use only

1.3. Details of the supplier of the safety data sheet

Veltek Associates, Inc. In Canada Distributed by:
15 Lee Blvd Canada Clean Room (CCR)
Malvern, PA 19355-1234 USA 20 Cope Dr.
T +1 610-644-8335 - F +1 610-644-8336 Kanata, ON K2M 2V8, Canada
vai@sterile.com Telephone: (888)595-8070

1.4. Emergency telephone number

Emergency number : CARECHEM 24 call: 1-866-928-0789

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification (GHS-US)

Not classified

2.2. Label elements

GHS-US labeling

No labeling applicable

2.3. Other hazards

No additional information available

2.4. Unknown acute toxicity (GHS-US)

No data available

SECTION 3: Composition/information on ingredients

3.1. Substance

Not applicable

Full text of H-phrases: see section 16

3.2. Mixture

Name	Product identifier	%	Classification (GHS-US)
Water	(CAS No) 7732-18-5	95-99.9	Not classified
Nitrogen	(CAS No) 7727-37-9	> 1	Not classified
Potassium tall oil soap	(CAS No) 61790-44-1	< 0.1	Eye Irrit. 2A, H319
Nonionic detergent		< 0.1	Not classified

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general : Never give anything by mouth to an unconscious person. If you feel unwell, seek medical advice (show the label where possible).
First-aid measures after inhalation : Remove victim to fresh air and keep at rest in a position comfortable for breathing. If symptoms develop obtain medical attention.
First-aid measures after skin contact : Rinse skin with water/shower. Remove/Take off immediately all contaminated clothing. Gently wash with plenty of soap and water. If symptoms develop, obtain medical attention.
First-aid measures after eye contact : Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If symptoms develop, obtain medical attention.

DECON-CLEAN[®] WFI FORMULA AEROSOL

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

First-aid measures after ingestion : Get immediate medical attention. Rinse mouth. Do NOT induce vomiting. Do not give an unconscious person anything to drink.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/injuries after skin contact : Repeated and/or prolonged skin contact may cause irritation.
Symptoms/injuries after eye contact : May cause eye irritation.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Dry powder. Carbon dioxide. Alcohol resistant foam.
Unsuitable extinguishing media : Do not use a heavy water stream.

5.2. Special hazards arising from the substance or mixture

Fire hazard : Not flammable.
Explosion hazard : Containers may rupture when heated.
Reactivity : Stable under normal conditions.

5.3. Advice for firefighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Exercise caution when fighting any chemical fire. Prevent fire-fighting water from entering environment.
Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection. Use self-contained breathing apparatus when in close proximity to fire.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Remove ignition sources. Ensure adequate ventilation.

6.1.1. For non-emergency personnel

Emergency procedures : Evacuate unnecessary personnel.

6.1.2. For emergency responders

Protective equipment : Equip cleanup crew with proper protection.
Emergency procedures : Ensure adequate ventilation. Avoid inhalation of vapors.

6.2. Environmental precautions

Prevent entry to sewers and public waters. Notify authorities if large amounts of the product enters sewers or public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Soak up spills with inert solids, such as clay or diatomaceous earth as soon as possible. Collect spillage. Store away from other materials.

6.4. Reference to other sections

SECTION 8: Exposure controls/personal protection. SECTION 13: Disposal considerations.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Provide good ventilation in process area to prevent formation of vapor. Containers may rupture when heated.
Hygiene measures : Do not eat, drink or smoke when using this product. Wash hands thoroughly after handling.

7.2. Conditions for safe storage, including any incompatibilities

Technical measures : Keep away from open flames, hot surfaces and sources of ignition.
Storage conditions : Keep in a cool, well-ventilated place away from incompatible materials.
Incompatible products : Strong bases. Strong acids. Direct sunlight. Heat sources.

7.3. Specific end use(s)

Residue Remover.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

DECON-CLEAN[®] WFI FORMULA AEROSOL

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

8.2. Exposure controls

Appropriate engineering controls	: Ensure adequate ventilation. Emergency eye wash fountains and safety showers should be available in the immediate vicinity of any potential exposure.
Personal protective equipment	: Avoid all unnecessary exposure. Wash hands thoroughly after handling.
Hand protection	: Wear suitable gloves.
Eye protection	: Chemical goggles or safety glasses.
Skin and body protection	: Wear suitable protective clothing.
Respiratory protection	: In case of insufficient ventilation, wear suitable respiratory equipment.
Thermal hazard protection	: Not required for normal conditions of use.
Environmental exposure controls	: Contain any spills with dikes or absorbents to prevent migration and entry into sewers or streams.
Other information	: Do not eat, drink or smoke during use. Do not breathe Vapors. Do not eat, drink or smoke when using this product. Spray away from eyes and face.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	: Liquid
Color	: Pale yellow.
Odor	: Slight.
Odor threshold	: No data available
pH	: 8 - 10.5
Relative evaporation rate (butyl acetate=1)	: No data available
Melting point	: No data available
Freezing point	: 32 °F (0 °C)
Boiling point	: 212 °F (100 °C)
Flash point	: Not applicable
Auto-ignition temperature	: Not applicable
Decomposition temperature	: No data available
Flammability (solid, gas)	: Not applicable
Vapor pressure	: Not applicable
Relative vapor density at 20 °C	: Not applicable
Relative density	: No data available
Solubility	: Miscible. Water: 100 %
Log Pow	: No data available
Log Kow	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: Not explosive.
Oxidizing properties	: Not oxidizing.
Explosive limits	: Not applicable

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under normal conditions.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

Containers may rupture when heated.

10.4. Conditions to avoid

Direct sunlight. Extremely high or low temperatures.

DECON-CLEAN[®] WFI FORMULA AEROSOL

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

10.5. Incompatible materials

Strong acids. Strong bases. Oxidizing agents.

10.6. Hazardous decomposition products

Carbon monoxide. Carbon dioxide.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity	:	Not classified
Skin corrosion/irritation	:	Not classified
		pH: 8 - 10.5
Serious eye damage/irritation	:	Not classified
		pH: 8 - 10.5
Respiratory or skin sensitization	:	Not classified
Germ cell mutagenicity	:	Not classified
		Based on available data, the classification criteria are not met
Carcinogenicity	:	Not classified
Reproductive toxicity	:	Not classified
		Based on available data, the classification criteria are not met
Specific target organ toxicity (single exposure)	:	Not classified
Specific target organ toxicity (repeated exposure)	:	Not classified
		Based on available data, the classification criteria are not met
Aspiration hazard	:	Not classified
		Based on available data, the classification criteria are not met
Potential Adverse human health effects and symptoms	:	Harmful if swallowed.
Symptoms/injuries after skin contact	:	Repeated and/or prolonged skin contact may cause irritation.
Symptoms/injuries after eye contact	:	May cause eye irritation.

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

DECON-CLEAN[®] WFI FORMULA AEROSOL	
Persistence and degradability	No data available.

12.3. Bioaccumulative potential

DECON-CLEAN[®] WFI FORMULA AEROSOL	
Bioaccumulative potential	Low bioaccumulation potential.

12.4. Mobility in soil

DECON-CLEAN[®] WFI FORMULA AEROSOL	
Ecology - soil	No data available.

12.5. Other adverse effects

Other information : Avoid release to the environment.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste disposal recommendations	:	Dispose in a safe manner in accordance with local/national regulations. Dispose of this material and its container to hazardous or special waste collection point.
Additional information	:	Handle empty containers with care.
Ecology - waste materials	:	Avoid release to the environment.

DECON-CLEAN[®] WFI FORMULA AEROSOL

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

SECTION 14: Transport information

In accordance with DOT

UN-No.(DOT) : 1950
Proper Shipping Name (DOT) : Aerosols
Hazard labels (DOT) : 2.2 - Non-flammable gas



DOT Packaging Exceptions (49 CFR 173.xxx) : 306
DOT Packaging Non Bulk (49 CFR 173.xxx) : None
DOT Packaging Bulk (49 CFR 173.xxx) : None

Additional information

Other information : Not classified.
Special transport precautions : No special precautions required.

Transport by sea

UN-No. (IMDG) : 1950
Proper Shipping Name (IMDG) : AEROSOLS
Class (IMDG) : 2 - Gases

Air transport

UN-No. (IATA) : 1950
Proper Shipping Name (IATA) : AEROSOLS
Class (IATA) : 2

SECTION 15: Regulatory information

15.1. US Federal regulations

Nitrogen (7727-37-9)

Listed on the United States TSCA (Toxic Substances Control Act) inventory

Water (7732-18-5)

Listed on the United States TSCA (Toxic Substances Control Act) inventory

Potassium tall oil soap (61790-44-1)

Listed on the United States TSCA (Toxic Substances Control Act) inventory

15.2. International regulations

CANADA

DECON-CLEAN[®] WFI FORMULA AEROSOL

WHMIS Classification : Uncontrolled product according to WHMIS classification criteria

15.2.2. National regulations

No additional information available

15.3. US State regulations

Nitrogen (7727-37-9)

U.S. - New Jersey - Right to Know Hazardous Substance List

DECON-CLEAN[®] WFI FORMULA AEROSOL

Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

SECTION 16: Other information

Data sources : US OSHA HazCom (GHS) 25 May 2012.


Full text of H-phrases: see section 16:

Eye Irrit. 2A	Serious eye damage/eye irritation Category 2A
H319	Causes serious eye irritation

NFPA health hazard : 0 - Exposure under fire conditions would offer no hazard beyond that of ordinary combustible materials.

NFPA fire hazard : 0 - Materials that will not burn.

NFPA reactivity : 1 - Normally stable, but can become unstable at elevated temperatures and pressures or may react with water with some release of energy, but not violently.



HMIS III Rating

Health : 0 Minimal Hazard - No significant risk to health

Flammability : 0 Minimal Hazard

Physical : 1 Slight Hazard

NCEC SDS US GHS (Hazcom 2012)

Information contained in this publication or as otherwise supplied to Users is believed to be accurate and is given in good faith, but it is for the Users to satisfy themselves of the suitability of the product for their own particular purpose. Veltek Associates, Inc. gives no warranty as to the fitness of the product for any particular purpose and any implied warranty or condition (statutory or otherwise) is excluded except to the extent that exclusion is prevented by law. Veltek Associates, Inc. accepts no liability for loss or damage resulting from reliance on this information. Freedom under Patents, Copyright and Designs cannot be assumed.