



# HYPO-CHLOR® 0.52% WIPE

## Safety Data Sheet

VELTEK ASSOCIATES, INC.

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations  
Issue date: 4/24/2015 Revision date: 9/16/2021 Supersedes: 5/22/2018 Version: 3.0

### SECTION 1: Identification

#### 1.1. Identification

Product form : Mixture  
Product name : HYPO-CHLOR® 0.52% WIPE  
Product code : SDS VEL-130

#### 1.2. Recommended use and restrictions on use

Use of the substance/mixture : Cleaning wipe  
For professional use only

#### 1.3. Supplier

Veltek Associates, Inc.  
15 Lee Blvd  
Malvern, PA 19355-1234 USA  
Telephone: +1 610-644-8335 - Fax: +1 610-644-8336  
E-mail: vai@sterile.com

In Canada distributed by:  
Canada Clean Room (CCR)  
20 Cope Dr.  
Kanata, ON K2M 2V8, Canada  
Telephone: (888)595-8070

#### 1.4. Emergency telephone number

Emergency number : CARECHEM 24: 1-215-207-0061  
1-866-928-0789 (toll free)  
Canada: 1-800-579-7421 (toll free)  
Mexico: +52-55-5004-8763

### SECTION 2: Hazard(s) identification

#### 2.1. Classification of the substance or mixture

##### GHS US classification

Corrosive to metals Category 1	H290	May be corrosive to metals
Hazardous to the aquatic environment - Acute Hazard Category 2	H401	Toxic to aquatic life
Hazardous to the aquatic environment - Chronic Hazard Category 3	H412	Harmful to aquatic life with long lasting effects

Full text of H statements : see section 16

#### 2.2. GHS Label elements, including precautionary statements

##### GHS US labeling

Hazard pictograms (GHS US) :



Signal word (GHS US) : Warning  
Hazard statements (GHS US) : H290 - May be corrosive to metals  
H401 - Toxic to aquatic life  
H412 - Harmful to aquatic life with long lasting effects

# HYPO-CHLOR<sup>®</sup> 0.52% WIPE

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Precautionary statements (GHS US) : P234 - Keep only in original container.  
P273 - Avoid release to the environment.  
P390 - Absorb spillage to prevent material-damage.  
P406 - Store in corrosive resistant container with a resistant inner liner.  
P501 - Dispose of contents/container to an authorized waste collection point.

### 2.3. Other hazards which do not result in classification

No additional information available

### 2.4. Unknown acute toxicity (GHS US)

Not applicable

## SECTION 3: Composition/Information on ingredients

### 3.1. Substances

Not applicable

### 3.2. Mixtures

Name	Product identifier	%	GHS US classification
Sodium hypochlorite	CAS-No.: 7681-52-9	0.40 - 0.65	Met. Corr. 1, H290 Skin Corr. 1B, H314 Eye Dam. 1, H318 STOT SE 3, H335 Aquatic Acute 1, H400 (M=10) Aquatic Chronic 1, H410

Full text of hazard classes and H-statements : see section 16

## SECTION 4: First-aid measures

### 4.1. Description of first aid measures

First-aid measures general : Never give anything by mouth to an unconscious person. If you feel unwell, seek medical advice (show the label where possible).  
First-aid measures after inhalation : Remove person to fresh air and keep at rest in a position comfortable for breathing. If symptoms develop obtain medical attention.  
First-aid measures after skin contact : Remove affected clothing and wash all exposed skin area with mild soap and water, followed by warm water rinse. If skin irritation occurs: Get medical advice/attention.  
First-aid measures after eye contact : Rinse immediately with plenty of water. Remove contact lenses, if present and easy to do. Continue rinsing. If symptoms develop, obtain medical attention.  
First-aid measures after ingestion : Do NOT induce vomiting. Rinse mouth. Never give anything by mouth to an unconscious person. If symptoms develop, obtain medical attention.

### 4.2. Most important symptoms and effects (acute and delayed)

Symptoms/effects after skin contact : Repeated and/or prolonged skin contact may cause irritation.  
Symptoms/effects after eye contact : Slight eye irritant upon direct contact.

### 4.3. Immediate medical attention and special treatment, if necessary

Treat symptomatically.

# HYPO-CHLOR<sup>®</sup> 0.52% WIPE

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

### SECTION 5: Fire-fighting measures

#### 5.1. Suitable (and unsuitable) extinguishing media

Suitable extinguishing media : Use extinguishing media appropriate for surrounding fire.

Unsuitable extinguishing media : None known.

#### 5.2. Specific hazards arising from the chemical

Fire hazard : Not flammable.

Hazardous decomposition products in case of fire : Carbon monoxide. Carbon dioxide.

#### 5.3. Special protective equipment and precautions for fire-fighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Exercise caution when fighting any chemical fire. Prevent fire-fighting water from entering environment.

Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection. Use self-contained breathing apparatus when in close proximity to fire.

### SECTION 6: Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

##### 6.1.1. For non-emergency personnel

Emergency procedures : Evacuate unnecessary personnel. Ventilate area. Avoid inhalation of vapors. Avoid contact with skin, eyes and clothing.

##### 6.1.2. For emergency responders

Protective equipment : Use personal protective equipment as required. See Section 8.

Emergency procedures : Ventilate area. Avoid inhalation of vapors. Avoid contact with eyes, skin and clothing.

#### 6.2. Environmental precautions

Do not allow to enter drains or water courses. Notify authorities if product enters sewers or public waters.

#### 6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Place in a suitable container for disposal in accordance with the waste regulations (see Section 13). Store in corrosive resistant container with a resistant inner liner. Store away from other materials.

#### 6.4. Reference to other sections

SECTION 8: Exposure controls/personal protection. SECTION 13: Disposal considerations.

### SECTION 7: Handling and storage

#### 7.1. Precautions for safe handling

Precautions for safe handling : Provide good ventilation in process area to prevent formation of vapor. Wear suitable protective clothing, gloves and eye or face protection. Avoid contact with skin, eyes and clothing. Avoid inhalation of vapors.

Hygiene measures : Handle in accordance with good industrial hygiene and safety practice. Do not eat, drink or smoke when using this product. Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Wash contaminated clothing before reuse.

#### 7.2. Conditions for safe storage, including any incompatibilities

Technical measures : Comply with applicable regulations.

# HYPO-CHLOR<sup>®</sup> 0.52% WIPE

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Storage conditions	: Keep only in the original container in a cool, well ventilated place away from : acids. Store in corrosive resistant container with a resistant inner liner. Keep container closed when not in use.
Incompatible materials	: Acids. Water-reactive materials. Steel.

### SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

No additional information available

#### 8.2. Appropriate engineering controls

Appropriate engineering controls	: Provide good ventilation in process area to prevent formation of vapor.
Environmental exposure controls	: Contain any spills with dikes or absorbents to prevent migration and entry into sewers or streams.

#### 8.3. Individual protection measures/Personal protective equipment

##### Personal protective equipment:

Avoid all unnecessary exposure.

<b>Hand protection:</b>
Wear chemically resistant protective gloves. The exact breakthrough time has to be found out by the manufacturer of the protective gloves and has to be observed. Gloves should be removed and replaced if there are any signs of degradation or breakthrough.
<b>Eye protection:</b>
Wear goggles or safety glasses with side shields if contact with the eyes is possible
<b>Skin and body protection:</b>
Long-sleeved protective clothing
<b>Respiratory protection:</b>
Not required for normal conditions of use. In case of insufficient ventilation, wear suitable respiratory equipment

##### Thermal hazard protection:

Not required for normal conditions of use.

##### Other information:

Do not eat, drink or smoke during use. Handle in accordance with good industrial hygiene and safety procedures.

### SECTION 9: Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

Physical state	: Liquid
Appearance	: Impregnated wipe.
Color	: Colorless to straw-colored liquid, white fabric
Odor	: Chlorine-like
Odor threshold	: No data available
pH	: 10.5 – 12.4 (Liquid)
Melting point	: Not applicable
Freezing point	: No data available
Boiling point	: 212 °F (100 °C)(Liquid)
Flash point	: Not applicable
Relative evaporation rate (butyl acetate=1)	: No data available
Flammability (solid, gas)	: Not flammable.
Vapor pressure	: No data available

# HYPO-CHLOR<sup>®</sup> 0.52% WIPE

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Relative vapor density at 20 °C	: No data available
Relative density	: 1 – 1.07 (Liquid) (Water = 1)
Solubility	: Water: Liquid: Miscible. Wipes: Insoluble
Log Pow	: No data available
Auto-ignition temperature	: Not applicable
Decomposition temperature	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosion limits	: No data available
Explosive properties	: Not explosive.
Oxidizing properties	: Not oxidizing.

### 9.2. Other information

No additional information available

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

Stable under recommended handling and storage conditions (see section 7).

### 10.2. Chemical stability

Slowly decomposes on contact with air.

### 10.3. Possibility of hazardous reactions

May produce small amounts of chlorine gas if mixed with incompatible materials.

### 10.4. Conditions to avoid

Extremely high or low temperatures.

### 10.5. Incompatible materials

Acids. Water-reactive materials. Steel.

### 10.6. Hazardous decomposition products

Carbon dioxide. Carbon monoxide. May produce small amounts of chlorine gas if mixed with incompatible materials.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity (oral)	: Not classified
Acute toxicity (dermal)	: Not classified
Acute toxicity (inhalation)	: Not classified

#### Sodium hypochlorite (7681-52-9)

LD50 oral, rat	8830 mg/kg (12.5% Aqueous solution)
LD50 dermal, rabbit	> 20000 mg/kg (12.5% Aqueous solution)

Skin corrosion/irritation	: Not classified pH: 10.5 – 12.4 (Liquid)
Serious eye damage/irritation	: Not classified pH: 10.5 – 12.4 (Liquid)
Respiratory or skin sensitization	: Not classified
Germ cell mutagenicity	: Not classified

# HYPO-CHLOR® 0.52% WIPE

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Carcinogenicity : Not classified  
Reproductive toxicity : Not classified  
STOT-single exposure : Not classified

Sodium hypochlorite (7681-52-9)	
STOT-single exposure	May cause respiratory irritation.
STOT-repeated exposure	: Not classified
Aspiration hazard	: Not classified
Viscosity, kinematic	: No data available
Symptoms/effects after skin contact	: Repeated and/or prolonged skin contact may cause irritation.
Symptoms/effects after eye contact	: Slight eye irritant upon direct contact.

## SECTION 12: Ecological information

### 12.1. Toxicity

Ecology - general : Toxic to aquatic life. Harmful to aquatic life with long lasting effects.

Sodium hypochlorite (7681-52-9)	
LC50 fish	0.06 mg/l - 96 Hours (freshwater fish)
EC50 Daphnia	0.141 mg/l - 48 Hours (Daphnia magna)
LC50 fish 2	0.032 mg/l - 96 Hours (marine water fish)
EC50 - Crustacea [2]	35 µg/l - 48 Hours (Ceriodaphnia dubia)
ErC50 algae	0.0499 mg/l - 72 Hours (Freshwater)
NOEC chronic fish	0.04 mg/l - 28 days (Menidia peninsulae)
NOEC chronic crustacea	0.007 mg/l - 15 days (estimated)

### 12.2. Persistence and degradability

Sodium hypochlorite (7681-52-9)	
Persistence and degradability	Not relevant for inorganic substances.

### 12.3. Bioaccumulative potential

Sodium hypochlorite (7681-52-9)	
Log Pow	-3.42 (20 °C, pH 12.5, Quantitative structure-activity relationship (QSAR))

### 12.4. Mobility in soil

Sodium hypochlorite (7681-52-9)	
Ecology - soil	Miscible with water.

### 12.5. Other adverse effects

Other information : Avoid release to the environment.

## SECTION 13: Disposal considerations

### 13.1. Disposal methods

Waste disposal recommendations : Dispose in a safe manner in accordance with local/national regulations.

# HYPO-CHLOR<sup>®</sup> 0.52% WIPE

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Ecology - waste materials : Avoid release to the environment.

### SECTION 14: Transport information

In accordance with Department of Transport / Transportation of Dangerous Goods / IMDG / IATA

#### 14.1. UN number

DOT NA No : UN3244  
UN-No. (TDG) : UN3244  
UN-No. (IMDG) : 3244  
UN-No. (IATA) : 3244

#### 14.2. UN proper shipping name

Proper Shipping Name (DOT) : Solids containing corrosive liquid, n.o.s. (Sodium hypochlorite)  
Proper Shipping Name (TDG) : SOLIDS CONTAINING CORROSIVE LIQUID, N.O.S. (Sodium hypochlorite)  
Proper Shipping Name (IMDG) : SOLIDS CONTAINING CORROSIVE LIQUID, N.O.S. (Sodium hypochlorite)  
Proper Shipping Name (IATA) : Solids containing corrosive liquid, n.o.s. (Sodium hypochlorite)  
Transport document description (DOT) : UN3244 Solids containing corrosive liquid, n.o.s. (Sodium hypochlorite), 8, II  
Transport document description (TDG) : UN3244 SOLIDS CONTAINING CORROSIVE LIQUID, N.O.S. (Sodium hypochlorite), 8, II  
Transport document description (IMDG) : UN 3244 SOLIDS CONTAINING CORROSIVE LIQUID, N.O.S. (Sodium hypochlorite), 8, II  
Transport document description (IATA) : UN 3244 Solids containing corrosive liquid, n.o.s. (Sodium hypochlorite), 8, II

#### 14.3. Transport hazard class(es)

##### DOT

Transport hazard class(es) (DOT) : 8  
Hazard labels (DOT) : 8



##### TDG

Transport hazard class(es) (TDG) : 8  
Hazard labels (TDG) : 8



##### IMDG

Transport hazard class(es) (IMDG) : 8  
Hazard labels (IMDG) : 8



##### IATA

Transport hazard class(es) (IATA) : 8  
Hazard labels (IATA) : 8

# HYPO-CHLOR® 0.52% WIPE

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations



### 14.4. Packing group

Packing group (DOT)	: II
Packing group (TDG)	: II
Packing group (IMDG)	: II
Packing group (IATA)	: II

### 14.5. Environmental hazards

Other information : No supplementary information available.

### 14.6. Special precautions for user

#### DOT

UN-No.(DOT)	: UN3244
DOT Special Provisions (49 CFR 172.102)	: 49, IB5, T3, TP33, TP3
DOT Packaging Exceptions (49 CFR 173.xxx)	: 154
DOT Packaging Non Bulk (49 CFR 173.xxx)	: 212
DOT Packaging Bulk (49 CFR 173.xxx)	: 240
DOT Quantity Limitations Passenger aircraft/rail (49 CFR 173.27)	: 15 kg
DOT Quantity Limitations Cargo aircraft only (49 CFR 175.75)	: 50 kg
DOT Vessel Stowage Location	: B
DOT Vessel Stowage Other	: 40

#### TDG

UN-No. (TDG)	: UN3244
TDG Special Provisions	: 16
Explosive Limit and Limited Quantity Index	: 1 kg
Excepted quantities (TDG)	: E2
Passenger Carrying Road Vehicle or Passenger Carrying Railway Vehicle Index	: 15 kg
Emergency Response Guide (ERG) Number	: 154

#### IMDG

No data available

#### IATA

No data available

### 14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

## SECTION 15: Regulatory information

### 15.1. US Federal regulations

#### HYPO-CHLOR® 0.52% WIPE

SARA Section 311/312 Hazard Classes	Physical hazard - Corrosive to metals
-------------------------------------	---------------------------------------



# HYPO-CHLOR<sup>®</sup> 0.52% WIPE

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

### Sodium hypochlorite (7681-52-9)

Listed on the United States TSCA (Toxic Substances Control Act) inventory

CERCLA RQ

100 lb

### 15.2. International regulations

#### CANADA

### Sodium hypochlorite (7681-52-9)

Listed on the Canadian DSL (Domestic Substances List)

#### EU-Regulations

No additional information available

#### National regulations

No additional information available

### 15.3. US State regulations

California Proposition 65 - This product does not contain any substances known to the state of California to cause cancer, developmental and/or reproductive harm

Component	State or local regulations
Sodium hypochlorite(7681-52-9)	U.S. - Massachusetts - Right To Know List; U.S. - New Jersey - Right to Know Hazardous Substance List; U.S. - Pennsylvania - RTK (Right to Know) List

## SECTION 16: Other information

Revision date : 09/16/2021  
Data sources : US OSHA HazCom (GHS) 25 May 2012.  
Other information : None.

### Full text of H-phrases

H290	May be corrosive to metals
H314	Causes severe skin burns and eye damage
H318	Causes serious eye damage
H335	May cause respiratory irritation
H400	Very toxic to aquatic life
H401	Toxic to aquatic life
H410	Very toxic to aquatic life with long lasting effects
H412	Harmful to aquatic life with long lasting effects

### Abbreviations and acronyms

	ACGIH (American Conference of Government Industrial Hygienists)
	CAS (Chemical Abstracts Service) number
	EC50 (Effective Concentration 50%)


# HYPO-CHLOR<sup>®</sup> 0.52% WIPE

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Abbreviations and acronyms	
	IARC (International Agency for Research on Cancer)
	IATA (International Air Transport Association)
	IMDG (International Maritime Dangerous Goods Code)
	IMO (International Maritime Organisation)
	LC50 (Lethal Concentration 50%)
	LD50 (Lethal Dose 50%)
	NOEC (No Observed Effect Concentration)
	OECD (Organisation for Economic Co-operation and Development)
	OSHA (Occupational Safety and Health Administration) (US)
	STEL (Short Term Exposure Limit)
	TSCA (Toxic Substances Control Act) (US)
	TWA (Time Weighted Average)
	UNxxxx (Number assigned by the United Nations Committee of Experts on the Transport of Dangerous Goods)

- NFPA health hazard : 1 - Materials that, under emergency conditions, can cause significant irritation.
- NFPA fire hazard : 0 - Materials that will not burn under typical fire conditions, including intrinsically noncombustible materials such as concrete, stone, and sand.
- NFPA reactivity : 1 - Materials that in themselves are normally stable but can become unstable at elevated temperatures and pressures.



- Hazard Rating
- Health : 1 Slight Hazard - Irritation or minor reversible injury possible
- Flammability : 0 Minimal Hazard - Materials that will not burn
- Physical : 1 Slight Hazard - Materials that are normally stable but can become unstable (self-react) at high temperatures and pressures. Materials may react non-violently with water or undergo hazardous polymerization in the absence of inhibitors.

Indication of changes:			
Section	Changed item	Change	Comments
2	Hazards identification	Modified	
5	Fire fighting measures	Modified	
5	Fire fighting measures	Modified	
6	Accidental release measures	Modified	
7	Handling and storage	Modified	
9	Physical and chemical properties	Modified	
10	Stability and reactivity	Modified	
11	Toxicological information	Modified	
12.	Ecological information	Modified	
14	Transport information	Modified	

# HYPO-CHLOR<sup>®</sup> 0.52% WIPE

## Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

15	Regulatory information	Modified	
16	Other information	Modified	

Safety Data Sheet (SDS), USA

Information contained in this publication or as otherwise supplied to Users is believed to be accurate and is given in good faith, but it is for the Users to satisfy themselves of the suitability of the product for their own particular purpose. Veltek Associates, Inc. gives no warranty as to the fitness of the product for any particular purpose and any implied warranty or condition (statutory or otherwise) is excluded except to the extent that exclusion is prevented by law. Veltek Associates, Inc. accepts no liability for loss or damage resulting from reliance on this information. Freedom under Patents, Copyright and Designs cannot be assumed.

This SDS has been translated into the official language of the country/region in which the product is to be placed on the market. Where no official translation exists, the regulatory text is reported in English, as it appears in the relevant regulatory text.